

Stomacare
product guide

unoro

2026

Here for what matters



Welcome to Unora

We're here to redefine what's possible in stoma care – for people living with change and for you, who share our commitment to helping them live the lives they deserve.

Our holistic and flexible portfolio of products and accessories is designed with protection and comfort in mind.

Combining innovative materials and decades of experience, our five material options make it easy to shape care to each person – whatever their needs and wherever they are on their journey.

And because care goes beyond the product itself, we're here for the quiet victories, the hard moments and everything in between. From the initial conversations to life after surgery, we're the partner and life-line, ready with guidance and support that helps everyone feel confident every step of the way.

With Unora, you can focus on what matters most – delivering the best care to more people.

We're here for what matters.

Contents

Introduction

Product families	8
Skin protection materials	10
Convexities	12
Systems	14
Product features	16

Products

Colostomy	18
Ileostomy	34
Urostomy	50
Special supply	66
Accessories	72

Table of contents by article number	78
-------------------------------------	----



Product families

Our portfolio spans solutions designed to support different stoma care needs. From skin protection products made from innovative silicone to traditional hydrocolloid materials – all products have been developed to ensure comfortable and reliable care in everyday life.



Silkura®

Breakthrough silicone technology for skin health.

Breathable silicone skin protection that moves naturally with the body, delivering a secure fit that feels light and comfortable throughout the day.

- Breathable silicone barrier helps maintain the skin's natural moisture balance
- Non-reactive, skin-friendly material suitable for sensitive or fragile skin
- Lightweight, water-resistant construction that remains dry against the skin
- Gentle, consistent adhesion designed to protect skin integrity over time
- Comfortable removal without damaging the skin
- Pouch with discreet, elegant black design, 100% opaque for confidence and discretion

Stomocur®

Reliable hydrocolloid technology for everyday confidence.

Stomocur® systems combine trusted adhesion with a range of hydrocolloid based skin protection materials that help adapt care to different skin and stoma needs.

- A selection of materials for every need, e.g., infused with alginate or materials free of animal-derived ingredients
- Proven adhesion for a secure and comfortable seal
- Designed to support peristomal skin integrity and long wear time
- Wide range of pouch shapes, sizes and convex options to support individual anatomy
- Available as a one-piece or in a two-piece click system

Skin protection materials

Protecting peristomal skin is central to effective stoma care. Our portfolio includes a range of skin protection materials designed to support different skin needs and wearing preferences, making it easier to share care around the individual.



Silicone

Breathable silicone skin protection that combines gentle adhesion with a secure seal.

- Breathable barrier that helps maintain the skin's natural moisture balance
- Adapts closely to the skin's surface to help create a secure seal
- Moves naturally with the body for a light, barely noticeable feel



Hydrocolloid

Trusted hydrocolloid skin protection with reliable adhesion and a secure seal to fit everyday needs.

- Established hydrocolloid formulation for dependable adhesion
- Absorbs moisture around the stoma to help maintain a secure seal
- Proven material technology widely used in everyday ostomy care



Combiflange

Dual-layer hydrocolloid technology engineered to enhance adhesion and security while remaining flexible and comfortable.

- Two hydrocolloid layers bonded together to increase absorptive capacity
- 0.4 mm thicker structure compared to standard hydrocolloid for enhanced adhesion
- Thin, tapered border designed to adapt to uneven skin contours



Alginate

Alginate-infused skin protection designed to absorb up to 1.5× more fluid than standard hydrocolloids for improved moisture management.

- Contains 15% alginate derived from natural seaweed
- Forms an elastic gel that adapts to any skin surface
- Helps maintain a balanced peristomal environment and supports comfort and extended wear time



Microskin

Ultra-thin, flexible skin protection that adapts closely to the body and supports natural movement.

- Ultra-thin foil with a thickness of just 0.02 mm for improved discretion
- Highly flexible design that adapts to body contours and movement
- Suitable for active lifestyles



Gelatine free


Hydrocolloid technology without gelatine or animal-derived components to support individual preferences.

- Reliable adhesion designed to maintain a secure seal
- Suitable for individuals with specific dietary, cultural or ethical preferences
- Comfortable skin protection that supports everyday wear

Convexities


Every stoma and body contour is different. When additional support is needed to help the stoma protrude and maintain a secure seal, convex baseplates can help create a more reliable fit.

Our convex portfolio offers different convexity heights and plateau sizes, making it easier to match the right level of support to each person's needs.



Medium convex

Medium convex baseplates provide gentle support while maintaining flexibility and comfort. With a convex height of 3.5 mm, they help create a more secure fit and adapt naturally to the body's contours. Medium convex baseplates are available in three different plateau sizes.



Convex

Convex baseplates provide a deeper level of support when additional assistance is needed to help the stoma protrude and maintain a secure seal. With a convex height of 7 mm, they are designed to support more challenging stoma profiles while maintaining comfort and flexibility. Convex baseplates are available in up to four different plateau sizes.

Available in Combiflange and Alginate materials, our Medium convex and Convex options combine reliable adhesion with soft, flexible construction designed for comfortable everyday wear.


Systems

Different pouching systems support different approaches to stoma care. Our one-piece and two-piece systems allow care to be adapted to individual preferences, skin needs and daily routines.



One-piece systems

In a one-piece system, the pouch and baseplate form a single unit. This creates a simple and discreet solution where the entire system is changed together.



Two-piece clic system

Two-piece systems allow the pouch to be changed without removing the baseplate from the skin. This can help reduce the frequency of barrier removal while offering greater flexibility in daily care.

The Unora two-piece system combines dependable adhesion with a secure locking ring connection, allowing pouches to be attached and removed with confidence.

- Easy-to-handle coupling system
- Audible "click" to confirm pouch attachment
- Locking ring available in three sizes (45 mm, 57 mm and 70 mm)


Product features

Our pouches are developed and manufactured in Berlin using modern production technologies and high-quality materials.


They are characterised by innovative and unique properties that meet the medical requirements of healthcare professionals and users.

Pouch shapes


Symmetric



Asymmetric




Outlets




Cover material

Non-woven cover



Textile cover




Filter

- Colostomy – Activated carbon filter
- Ileostomy – MPC filter
MPC filter technology (Moisture Protect Composition)
= additional protective layer that prevents liquid absorption.

Window

All of our pouches feature a window invented by Unora. This enables easy positioning and observation of the stoma.



Colostomy

Unora colostomy solutions combine reliable filters, comfortable materials and skin-friendly barriers to help protect the peristomal area and maintain confidence in daily wear.

The Unora colostomy portfolio spans two complementary product families that make it easier to shape care around individual anatomy, skin condition and lifestyle – whether supporting established routines or helping someone regain confidence after surgery.

Silkura®

Breakthrough silicone technology for skin health.



Silicone

Breathable silicone skin protection that combines gentle adhesion with a secure seal.

Stomocur®

Reliable hydrocolloid technology for everyday confidence.



Hydrocolloid

Trusted hydrocolloid skin protection with reliable adhesion and a secure seal to fit everyday needs.



Combiflange

Dual-layer hydrocolloid technology engineered to enhance adhesion and security while remaining flexible and comfortable.



Alginate

Alginate-infused skin protection designed to absorb up to 1.5× more fluid than standard hydrocolloids for improved moisture management.



Gelatine free

Hydrocolloid technology without gelatine or animal-derived components to support individual preferences.

One-piece systems

Silkura®	22
Stomocur® Select Large	22
Stomocur® Select	23
Stomocur® Select Small	23
Stomocur® Select Extra Small	24
Stomocur® Select Extra Small transparent	24
Stomocur® Select Cap	24
Stomocur® Select Alginat	25
Stomocur® Select Alginat Small	25
Stomocur® Select gelatinefrei	25
Stomocur® Primo medium konvex	26
Stomocur® Select medium konvex	26
Stomocur® Primo Alginat medium konvex	27
Stomocur® Primo konvex	27
Stomocur® Select konvex	28
Stomocur® Primo Alginat konvex	28
Stomocur® Select Alginat konvex	29

Two-piece systems

Stomocur® Clic	30
Stomocur® Clic Alginat	30
Stomocur® Clic konvex	31
Stomocur® Clic Alginat konvex	31
Stomocur® Select Clic	32



Silkura®

New for 2026

Window, filter, ultra-soft, opaque and water-repellent cover on both sides, adhesive surface max. 130 mm, pouch size 154 x 202 mm, capacity 330 ml, package contents 30 pieces.



Silicone

Stoma opening	Article number
15 – 60 mm	CSIL15B



Stomocur® Select Large

Can be cut up to 90 mm

Window, filter, soft cover both sides, adhesive surface max. 137 mm, pouch size 169 x 224 mm, capacity 470 ml, package contents 30 pieces.



Hydrocolloid

Stoma opening	Article number
20 – 90 mm	CSH20L
25 mm	CSH25L
30 mm	CSH30L
35 mm	CSH35L
40 mm	CSH40L



Stomocur® Select

Window, filter, soft cover both sides, adhesive surface max. 111 mm, pouch size 155 x 210 mm, capacity 380 ml, package contents 30 pieces.



Hydrocolloid

Stoma opening	Article number
20 – 60 mm	CSH20
25 mm	CSH25
30 mm	CSH30
35 mm	CSH35
40 mm	CSH40



Stomocur® Select Small

Window, filter, soft cover both sides, adhesive surface max. 111 mm, pouch size 140 x 194 mm, capacity 250 ml, package contents 30 pieces.



Hydrocolloid

Stoma opening	Article number
20 – 60 mm	CSH20S
25 mm	CSH25S
30 mm	CSH30S
35 mm	CSH35S
40 mm	CSH40S



Stomocur® Select Extra Small

Also suitable
for child care

Window, filter, soft cover both sides, adhesive surface 80 mm (Ø), pouch size 115 x 160 mm, capacity 140 ml, package contents 30 pieces.



Hydrocolloid

Stoma opening	Article number
10 – 50 mm	CSH10XS

Select Extra Small transparent


Transparent, filter, soft cover one side.

Stoma opening	Article number
10 – 50 mm	CSH10XST

Stomocur® Select Cap

After
irrigation

Package contents 30 pieces.



Stoma opening	Article number	Cap size	Adhesive surface Ø	Filter
20 – 40 mm	CSC20SNF	60 mm	50 mm	No
20 – 50 mm	CSC20	75 mm	60 mm	Yes
25 – 50 mm	CSC25NF	75 mm	60 mm	No
30 – 60 mm	CSC30LNF	85 mm	70 mm	No



Hydrocolloid



Stomocur® Select Alginat

Window, filter, soft cover both sides, adhesive surface max. 111 mm, pouch size 155 x 210 mm, capacity 380 ml, package contents 30 pieces.



Alginate

Stoma opening	Article number
20 – 60 mm	CSA20
25 mm	CSA25
30 mm	CSA30
35 mm	CSA35
40 mm	CSA40



Stomocur® Select Alginat Small

Window, filter, soft cover both sides, adhesive surface max. 111 mm, pouch size 140 x 194 mm, capacity 250 ml, package contents 30 pieces.



Alginate

Stoma opening	Article number
20 – 60 mm	CSA20S
25 mm	CSA25S
30 mm	CSA30S
35 mm	CSA35S
40 mm	CSA40S




Stomocur® Select gelatinefrei

Window, filter, soft cover both sides, adhesive surface max. 110 mm, pouch size 155 x 210 mm, capacity 380 ml, package contents 30 pieces.



Gelatine free

Stoma opening	Article number
20 – 60 mm	CSG20



Stomocur® Primo medium konvex

3.5 mm
plateau height

Window, filter, ultra-soft, opaque and water-repellent cover both sides, adhesive surface 116 mm (∅), pouch size 154 x 202 mm, capacity 330 ml, package contents 30 pieces.



Combiflange

Stoma opening	Plateau size	Article number
20 – 47 mm	55 mm	CAH20M1B
15 – 40 mm	45 mm	CAH15M2B
10 – 28 mm	35 mm	CAH10M3B



Stomocur® Select medium konvex

3.5 mm
plateau height

Window, filter, soft cover both sides, adhesive surface 118 mm (∅), pouch size 145 x 205 mm, capacity 380 ml, package contents 30 pieces.



Combiflange

Stoma opening	Plateau size	Article number
20 – 47 mm	55 mm	CSH20M1
25 mm		CSH25M1
30 mm		CSH30M1
35 mm		CSH35M1
40 mm	45 mm	CSH40M1
15 – 40 mm		CSH15M2
20 mm		CSH20M2
25 mm	CSH25M2	
10 – 28 mm	35 mm	CSH10M3



Stomocur® Primo Alginat medium konvex

3.5 mm
plateau height

Window, filter, ultra-soft, opaque and water-repellent cover both sides, adhesive surface 116 mm (∅), pouch size 154 x 202 mm, capacity 330 ml, package contents 30 pieces.



Alginat

Stoma opening	Plateau size	Article number
20 – 47 mm	55 mm	CAA20M1B
15 – 40 mm	45 mm	CAA15M2B
10 – 28 mm	35 mm	CAA10M3B



Stomocur® Primo konvex

7 mm
plateau height

Window, filter, ultra-soft, opaque and water-repellent cover both sides, adhesive surface 116 mm (∅), pouch size 154 x 202 mm, capacity 330 ml, package contents 30 pieces.



Combiflange

Stoma opening	Plateau size	Article number
20 – 43 mm	50 mm	CAH20V1B
15 – 35 mm	40 mm	CAH15V2B
10 – 25 mm	30 mm	CAH10V3B



Stomocur® Select konvex

7 mm
plateau height

Window, filter, soft cover both sides, adhesive surface 116 mm (∅), pouch size 155 x 210 mm, capacity 380 ml, package contents 30 pieces.



CombiFlange

Stoma opening	Plateau size	Article number
20 – 57 mm	60 mm	CSH20V0
20 – 43 mm	50 mm	CSH20V1
25 mm		CSH25V1
30 mm		CSH30V1
35 mm		CSH35V1
40 mm		CSH40V1
15 – 35 mm	40 mm	CSH15V2
20 mm		CSH20V2
25 mm		CSH25V2
30 mm		CSH30V2
10 – 25 mm	30 mm	CSH10V3
15 mm		CSH15V3
20 mm		CSH20V3



Stomocur® Select Alginat konvex

7 mm
plateau height

Window, filter, soft cover both sides, adhesive surface 114 mm (∅), pouch size 155 x 210 mm, capacity 380 ml, package contents 30 pieces.



Alginat

Stoma opening	Plateau size	Article number
20 – 43 mm	50 mm	CSA20V1
25 mm		CSA25V1
30 mm		CSA30V1
35 mm		CSA35V1
40 mm		CSA40V1
15 – 35 mm	40 mm	CSA15V2
20 mm		CSA20V2
25 mm		CSA25V2
30 mm		CSA30V2
10 – 25 mm	30 mm	CSA10V3
15 mm		CSA15V3
20 mm		CSA20V3



Stomocur® Primo Alginat konvex

7 mm
plateau height

Window, filter, ultra-soft, opaque and water-repellent cover both sides, adhesive surface 114 mm (∅), pouch size 154 x 202 mm, capacity 330 ml, package contents 30 pieces.



Alginat

Stoma opening	Plateau size	Article number
20 – 43 mm	50 mm	CAA20V1B
15 – 35 mm	40 mm	CAA15V2B
10 – 25 mm	30 mm	CAA10V3B



Stomocur®

Clic

Plate size 116 mm (Ø), package contents 5 pieces.



Locking ring	Stoma opening	Article number
45 mm	15 – 35 mm	CK4515
	20 mm	CK4520
	25 mm	CK4525
	30 mm	CK4530
57 mm	15 – 40 mm	CK5715
	35 mm	CK5735
	40 mm	CK5740
70 mm	15 – 50 mm	CK7015
	45 mm	CK7045
	50 mm	CK7050



Hydrocolloid



Stomocur®

Clic Alginat

Plate size 114 mm (Ø), package contents 5 pieces.



Locking ring	Stoma opening	Article number
45 mm	15 – 35 mm	CA4515
57 mm	15 – 40 mm	CA5715



Alginate



Stomocur®

Clic konvex

7 mm
plateau height

Plate size 116 mm (Ø), package contents 3 pieces.



Locking ring	Stoma opening	Plateau size	Article number
45 mm	15 – 35 mm	40 mm	CK4515V2
	20 mm		CK4520V2
	25 mm		CK4525V2
	30 mm		CK4530V2
	15 – 25 mm	30 mm	CK4515V3
57 mm	15 – 43 mm	50 mm	CK5715V1
	20 mm		CK5720V1
	25 mm		CK5725V1
	30 mm		CK5730V1
	35 mm		CK5735V1
	15 – 35 mm	40 mm	CK5715V2
	15 – 25 mm	30 mm	CK5715V3



Combiflange



Stomocur®

Clic Alginat konvex

7 mm
plateau height

Plate size 114 mm (Ø), package contents 3 pieces.



Locking ring	Stoma opening	Plateau size	Article number
45 mm	15 – 35 mm	40 mm	CA4515V2
	15 – 25 mm	30 mm	CA4515V3
57 mm	15 – 43 mm	50 mm	CA5715V1
	15 – 35 mm	40 mm	CA5715V2
	15 – 25 mm	30 mm	CA5715V3



Alginate



Stomocur®

Select Clic

Window, filter, soft cover both sides, pouch size 155 x 210 mm, capacity 380 ml, package contents 30 pieces.

Locking ring	Article number
45 mm	CC45
57 mm	CC57
70 mm	CC70



Ileostomy

Ileostomy care often places greater demands on both products and skin. Higher output, moisture and frequent emptying mean that reliable sealing, skin protection and ease of use are essential.

Our ileostomy solutions are designed to support throughout the day. Advanced filter technology, secure closures and skin-friendly materials work together to protect the peristomal skin and help reduce leakage risk.

Our ileostomy portfolio covers two complementary product families offering a range of materials and pouch options to help adapt care to different anatomies, skin needs and daily routines.

Silkura®

Breakthrough silicone technology for skin health.



Silicone

Breathable silicone skin protection that combines gentle adhesion with a secure seal.

Stomocur®

Reliable hydrocolloid technology for everyday confidence.



Hydrocolloid

Trusted hydrocolloid skin protection with reliable adhesion and a secure seal to fit everyday needs.



Combiflange

Dual-layer hydrocolloid technology engineered to enhance adhesion and security while remaining flexible and comfortable.



Alginate

Alginate-infused skin protection designed to absorb up to 1.5x more fluid than standard hydrocolloids for improved moisture management.



Gelatin free

Hydrocolloid technology without gelatin or animal-derived components to support individual preferences.

One-piece systems

Silkura®	38
Stomocur® Protect	38
Stomocur® Protect Small	39
Stomocur® Protect Alginat	39
Stomocur® Protect gelatinefrei	40
Stomocur® Drainage Protect	40
Stomocur® Primo medium konvex	41
Stomocur® Protect medium konvex	41
Stomocur® Protect medium konvex transparent	41
Stomocur® Primo Alginat medium konvex	42
Stomocur® Drainage Protect medium konvex	42
Stomocur® Primo konvex	43
Stomocur® Protect konvex	43
Stomocur® Primo Alginat konvex	44
Stomocur® Protect Alginat konvex	44
Stomocur® Drainage Protect konvex	45
Stomocur® Drainage Protect Bed Bag	45

Two-piece systems

Stomocur® Clic	46
Stomocur® Clic Alginat	46
Stomocur® Clic konvex	47
Stomocur® Clic Alginat konvex	47
Stomocur® Protect Clic	48
Stomocur® Protect Clic Small	48
Stomocur® Drainage Protect Clic transparent	49



Silkura®

New for 2026

Window, filter, ultra-soft, opaque and water-repellent cover on both sides, adhesive surface max. 130 mm, pouch size 155 x 215 mm (fully closed), capacity 460 ml, package contents 30 pieces.

Stoma opening	Article number
15 – 60 mm	ISIL15B



Silicone



Stomocur® Protect

Window, filter, soft cover both sides, "Safe and Easy" closure, adhesive surface max. 111 mm, pouch size 157 x 215 mm (fully closed), capacity 530 ml, package contents 30 pieces.

Stoma opening	Article number
20 – 70 mm	IPH2070
15 – 60 mm	IPH15
20 mm	IPH20
25 mm	IPH25
30 mm	IPH30
35 mm	IPH35
40 mm	IPH40



Hydrocolloid

Protect Transparent

Transparent, filter, soft cover one side, "Safe and Easy" closure.

Stoma opening	Article number
15 – 60 mm	IPH15T



Stomocur® Protect Small

Also suitable for child care

Window, filter, soft cover both sides, "Safe and Easy" closure, adhesive surface max. 111 mm, pouch size 157 x 185 mm (fully closed), capacity 360 ml, package contents 30 pieces.

Stoma opening	Article number
15 – 60 mm	IPH15S



Hydrocolloid



Stomocur® Protect Alginat

Window, filter, soft cover both sides, "Safe and Easy" closure, adhesive surface max. 111 mm, pouch size 157 x 215 mm (fully closed), capacity 530 ml, package contents 30 pieces.

Stoma opening	Article number
15 – 60 mm	IPA15
20 mm	IPA20
25 mm	IPA25
30 mm	IPA30
35 mm	IPA35
40 mm	IPA40



Alginate



Stomocur® Protect gelatinefrei

Window, filter, soft cover both sides, "Safe and Easy" closure, adhesive surface max. 110 mm, pouch size 157 x 215 mm (fully closed), capacity 530 ml, package contents 30 pieces.

Stoma opening	Article number
20 – 60 mm	IPG20



Gelatine free



Stomocur® Drainage Protect

Window, filter, soft cover both sides, large soft outlet, adhesive surface 116 mm (Ø), pouch size 157 x 293 mm, capacity 540 ml, package contents 10 pieces.

Stoma opening	Article number
20 – 60 mm	DP5020



Combiflange



Stomocur® Primo medium konvex

3.5 mm
plateau height

Window, filter, ultra-soft, opaque and water-repellent cover both sides, "Safe and Easy" closure, adhesive surface 116 mm (Ø), pouch size 155 x 215 mm (fully closed), capacity 460 ml, package contents 10 pieces.

Stoma opening	Plateau size	Article number
15 – 47 mm	55 mm	IAH15M1B
15 – 40 mm	45 mm	IAH15M2B
10 – 28 mm	35 mm	IAH10M3B



Combiflange



Stomocur® Protect medium konvex

3.5 mm
plateau height

Window, filter, soft cover both sides, "Safe and Easy" closure, adhesive surface 116 mm (Ø), pouch size 157 x 215 mm (fully closed), capacity 530 ml, package contents 10 pieces.

Stoma opening	Plateau size	Article number
15 – 47 mm	55 mm	IPH15M1
20 mm		IPH20M1
25 mm		IPH25M1
30 mm		IPH30M1
35 mm		IPH35M1
40 mm		IPH40M1
15 – 40 mm	45 mm	IPH15M2
20 mm		IPH20M2
25 mm		IPH25M2
10 – 28 mm	35 mm	IPH10M3



Combiflange

Protect medium konvex transparent

Transparent, filter, soft cover one side, "Safe and Easy" closure.

Stoma opening	Plateau size	Article number
15 – 47 mm	55 mm	IPH15M1T
15 – 40 mm	45 mm	IPH15M2T
10 – 28 mm	35 mm	IPH10M3T



Stomocur® Primo Alginat medium konvex

3.5 mm
plateau height

Window, filter, ultra-soft, opaque and water-repellent cover both sides, "Safe and Easy" closure, adhesive surface 114 mm (Ø), pouch size 155 x 215 mm (fully closed), capacity 460 ml, package contents 10 pieces.



Alginat

Stoma opening	Plateau size	Article number
15 – 47 mm	55 mm	IAA15M1B
15 – 40 mm	45 mm	IAA15M2B
10 – 28 mm	35 mm	IAA10M3B



Stomocur® Drainage Protect medium konvex

3.5 mm
plateau height

Window, filter, soft cover both sides, large soft outlet, adhesive surface 116 mm (Ø), pouch size 157 x 293 mm, capacity 540 ml, package contents 10 pieces.



Combiflange

Stoma opening	Plateau size	Article number
20 – 47 mm	55 mm	DP20M1
15 – 40 mm	45 mm	DP15M2
10 – 28 mm	35 mm	DP10M3



Stomocur® Primo konvex

7 mm
plateau height

Window, filter, ultra-soft, opaque and water-repellent cover both sides, "Safe and Easy" closure, adhesive surface 116 mm (Ø), pouch size 155 x 215 mm (fully closed), capacity 460 ml, package contents 10 pieces.



Combiflange

Stoma opening	Plateau size	Article number
15 – 43 mm	50 mm	IAH15V1B
15 – 35 mm	40 mm	IAH15V2B
10 – 25 mm	30 mm	IAH10V3B



Stomocur® Protect konvex

7 mm
plateau height

Window, filter, soft cover both sides, "Safe and Easy" closure, adhesive surface 116 mm (Ø), pouch size 157 x 215 mm (fully closed), capacity 530 ml, package contents 10 pieces.



Combiflange

Stoma opening	Plateau size	Article number
20 – 57 mm	60 mm	IPH20V0
15 – 43 mm	50 mm	IPH15V1
20 mm		IPH20V1
25 mm		IPH25V1
30 mm		IPH30V1
35 mm		IPH35V1
40 mm	IPH40V1	
15 – 35 mm	40 mm	IPH15V2
20 mm		IPH20V2
25 mm		IPH25V2
30 mm		IPH30V2
10 – 25 mm	30 mm	IPH10V3
15 mm		IPH15V3
20 mm		IPH20V3



Alginat

Stomocur® Primo Alginat konvex

7 mm
plateau height

Window, filter, ultra-soft, opaque and water-repellent cover both sides, "Safe and Easy" closure, adhesive surface 114 mm (∅), pouch size 155 x 215 mm (fully closed), capacity 460 ml, package contents 10 pieces.

Stoma opening	Plateau size	Article number
15 – 43 mm	50 mm	IAA15V1B
15 – 35 mm	40 mm	IAA15V2B
10 – 25 mm	30 mm	IAA10V3B



Alginat

Stomocur® Protect Alginat konvex

7 mm
plateau height

Window, filter, ultra-soft, opaque and water-repellent cover both sides, "Safe and Easy" closure, adhesive surface 114 mm (∅), pouch size 155 x 215 mm (fully closed), capacity 460 ml, package contents 10 pieces.

Stoma opening	Plateau size	Article number
15 – 43 mm	50 mm	IPA15V1
20 mm		IPA20V1
25 mm		IPA25V1
30 mm		IPA30V1
35 mm		IPA35V1
40 mm	IPA40V1	
15 – 35 mm	40 mm	IPA15V2
10 – 25 mm	30 mm	IPA10V3



Combiflange

Stomocur® Drainage Protect konvex

7 mm
plateau height

Window, filter, soft cover both sides, large soft outlet, adhesive surface 116 mm (∅), pouch size 157 x 293 mm, capacity 540 ml, package contents 10 pieces.

Stoma opening	Plateau size	Article number
20 – 43 mm	50 mm	DP20V1
15 – 35 mm	40 mm	DP15V2
10 – 25 mm	30 mm	DP10V3



Stomocur® Drainage Protect Bed Bag

Large-lumen and kink-resistant tube, integrated adapter for connection to all Unora Drainage Protect pouches, package contents 5 pieces.

Tube length	Capacity	Article number
130 mm	2000 ml	DPB2000

Compatible with:

- Stomocur® Drainage Protect
- Stomocur® Drainage Protect konvex
- Stomocur® Drainage Protect medium konvex
- Stomocur® Drainage Protect Clic transparent
- Stomocur® Drainage Extra "Access Port"



Stoma opening	Plateau size	Article number
20 – 60 mm	–	DP5020
20 – 43 mm	50 mm	DP20V1
15 – 35 mm	40 mm	DP15V2
10 – 25 mm	30 mm	DP10V3
20 – 47 mm	55 mm	DP20M1
15 – 40 mm	45 mm	DP15M2
10 – 28 mm	35 mm	DP10M3
Locking ring: 70 mm	–	CDP70T



Stomocur®

Clic

Plate size 116 mm (∅), package contents 5 pieces.



Locking ring	Stoma opening	Article number
45 mm	15 – 35 mm	CK4515
	20 mm	CK4520
	25 mm	CK4525
	30 mm	CK4530
57 mm	15 – 40 mm	CK5715
	35 mm	CK5735
	40 mm	CK5740
70 mm	15 – 50 mm	CK7015
	45 mm	CK7045
	50 mm	CK7050



Hydrocolloid



Stomocur®

Clic Alginat

Plate size 114 mm (∅), package contents 5 pieces.



Locking ring	Stoma opening	Article number
45 mm	15 – 35 mm	CA4515
57 mm	15 – 40 mm	CA5715



Alginate



Stomocur®

Clic konvex

7 mm
plateau height

Plate size 116 mm (∅), package contents 3 pieces.



Locking ring	Stoma opening	Plateau size	Article number
45 mm	15 – 35 mm	40 mm	CK4515V2
	20 mm		CK4520V2
	25 mm		CK4525V2
	30 mm		CK4530V2
	15 – 25 mm	30 mm	CK4515V3
57 mm	15 – 43 mm	50 mm	CK5715V1
	20 mm		CK5720V1
	25 mm		CK5725V1
	30 mm		CK5730V1
	35 mm		CK5735V1
	15 – 35 mm	40 mm	CK5715V2
	15 – 25 mm	30 mm	CK5715V3



Combiflange



Stomocur®

Clic Alginat konvex

7 mm
plateau height

Plate size 114 mm (∅), package contents 3 pieces.



Locking ring	Stoma opening	Plateau size	Article number
45 mm	15 – 35 mm	40 mm	CA4515V2
	15 – 25 mm	30 mm	CA4515V3
57 mm	15 – 43 mm	50 mm	CA5715V1
	15 – 35 mm	40 mm	CA5715V2
	15 – 25 mm	30 mm	CA5715V3



Alginate



Stomocur®

Protect Clic

Window, filter, soft cover both sides, "Safe and Easy" closure, pouch size 157 x 215 mm (fully closed), capacity 530 ml, package contents 30 pieces.

Locking ring	Article number
45 mm	CI45
57 mm	CI57
70 mm	CI70



Stomocur®

Protect Clic Small

Window, filter, soft cover both sides, "Safe and Easy" closure, pouch size 157 x 185 mm (fully closed), capacity 360 ml, package contents 30 pieces.

Locking ring	Article number
45 mm	CI45S



Stomocur®

Drainage Protect Clic transparent

Transparent, filter, soft cover one side, large soft outlet, pouch size 157 x 293 mm, capacity 540 ml, package contents 30 pieces.

Locking ring	Article number
70 mm	CDP70T

Urostomy

Continuous drainage requires solutions that provide dependable flow while protecting the skin around the stoma. Comfort, security and simplicity are key in everyday urostomy care.

Unora urostomy solutions combine reliable drainage with flexible skin barriers and adaptability for drainage systems, helping protect the skin and maintain confidence.

Our urostomy range is represented by the Stomocur® product family. All stoma products feature a non-twisting pouch design and a non-return valve. For the outlet, there is a flexible choice between a 'Safe' soft outlet or a 'Twist and Safe' tap valve for comfortable wear in everyday life.

Stomocur®

Reliable hydrocolloid technology for everyday confidence.



Combiflange

Dual-layer hydrocolloid technology engineered to enhance adhesion and security while remaining flexible and comfortable.



Alginate

Alginate-infused skin protection designed to absorb up to 1.5x more fluid than standard hydrocolloids for improved moisture management.



Microskin

Ultra-thin, flexible skin protection that adapts closely to the body and supports natural movement.



Gelatine free

Hydrocolloid technology without gelatine or animal-derived components to support individual preferences.

One-piece systems

Stomocur® Safe	54
Stomocur® Safe transparent	54
Stomocur® Twist & Safe	54
Stomocur® Twist & Safe transparent	54
Stomocur® Safe Microskin	55
Stomocur® Safe Microskin transparent	55
Stomocur® Twist & Safe Microskin	55
Stomocur® Twist & Safe Microskin transparent	55
Stomocur® Safe gelatinefrei transparent	56
Stomocur® Safe medium konvex	56
Stomocur® Safe medium konvex transparent	56
Stomocur® Twist & Safe medium konvex	57
Stomocur® Twist & Safe medium konvex transparent	57
Stomocur® Safe konvex	58
Stomocur® Safe konvex transparent	58
Stomocur® Twist & Safe konvex	59
Stomocur® Twist & Safe konvex transparent	59
Konticur® Urostomy Bed Bag	60
Stomocur® Uro Adapter	60
Stomocur® Drip Cap	60

Two-piece systems

Stomocur® Clic	62
Stomocur® Clic Alginat	62
Stomocur® Clic konvex	63
Stomocur® Clic Alginat konvex	63
Stomocur® Safe Clic	64
Stomocur® Safe Clic transparent	64
Stomocur® Twist & Safe Clic	64
Stomocur® Twist & Safe Clic transparent	64



Stomocur® Safe

Window, soft cover both sides, soft outlet, adhesive surface 116 mm (∅), pouch size 155 x 210 mm, capacity 330 ml, package contents 30 pieces with 5 adapters.

Stoma opening	Article number
10 – 60 mm	USH10F
15 mm	USH15F
20 mm	USH20F
25 mm	USH25F
30 mm	USH30F



Combiflange

Safe transparent

Transparent, soft cover one side, soft outlet.

Stoma opening	Article number
10 – 60 mm	USH10FT



Stomocur® Twist & Safe

Window, soft cover both sides, tap valve, adhesive surface 116 mm (∅), pouch size 155 x 210 mm, capacity 330 ml, package contents 30 pieces with 5 adapters & 10 drip caps.

Stoma opening	Article number
10 – 60 mm	USH10H
15 mm	USH15H
20 mm	USH20H
25 mm	USH25H
30 mm	USH30H



Combiflange

Twist & Safe transparent

Transparent, soft cover one side, tap valve.

Stoma opening	Article number
10 – 60 mm	USH10HT



Stomocur® Safe Microskin

Window, soft cover both sides, soft outlet, adhesive surface 101 mm x 114 mm, pouch size 155 x 210 mm, capacity 330 ml, package contents 30 pieces with 5 adapters.

Stoma opening	Article number
13 – 50 mm	USM13F



Microskin

Safe Microskin transparent

Transparent, soft cover one side, soft outlet.

Stoma opening	Article number
13 – 50 mm	USM13FT



Stomocur® Twist & Safe Microskin

Window, soft cover both sides, tap valve, adhesive surface 101 mm x 114 mm, pouch size 155 x 210 mm, capacity 330 ml, package contents 30 Piece with 5 adapters & 10 drip caps.

Stoma opening	Article number
13 – 50 mm	USM13H



Microskin

Twist & Safe Microskin transparent

Transparent, soft cover one side, tap valve.

Stoma opening	Article number
13 – 50 mm	USM13HT



Stomocur® Safe gelatinefrei transparent

Transparent, soft cover one side, soft outlet, adhesive surface max. 110 mm, pouch size 155 x 210 mm, capacity 330 ml, package contents 30 pieces with 5 adapters.

Stoma opening	Article number
15 – 60 mm	USG15FT



Gelatine free



Stomocur® Safe medium konvex

Window, soft cover both sides, soft outlet, adhesive surface 116 mm (Ø), pouch size 155 x 210 mm, capacity 330 ml, package contents 10 pieces with 2 adapters.

Stoma opening	Plateau size	Article number
10 – 47 mm	55 mm	USH10M1F
15 mm		USH15M1F
20 mm		USH20M1F
25 mm		USH25M1F
30 mm		USH30M1F
35 mm		USH35M1F
10 – 40 mm	45 mm	USH10M2F
15 mm		USH15M2F
20 mm		USH20M2F
25 mm		USH25M2F
10 – 28 mm	35 mm	USH10M3F



Combiflange

Safe medium konvex transparent

Transparent, soft cover one side, soft outlet.

Stoma opening	Plateau size	Article number
10 – 47 mm	55 mm	USH10M1FT
10 – 40 mm	45 mm	USH10M2FT
10 – 28 mm	35 mm	USH10M3FT



Stomocur® Twist & Safe medium konvex

3.5 mm
plateau height

Window, soft cover both sides, tap valve, adhesive surface 116 mm (Ø), pouch size 155 x 210 mm, capacity 330 ml, package contents 10 Set with 5 adapters & 5 drip caps.

Stoma opening	Plateau size	Article number
10 – 47 mm	55 mm	USH10M1H
15 mm		USH15M1H
20 mm		USH20M1H
25 mm		USH25M1H
30 mm		USH30M1H
35 mm		USH35M1H
10 – 40 mm	45 mm	USH10M2H
15 mm		USH15M2H
20 mm		USH20M2H
25 mm		USH25M2H
10 – 28 mm	35 mm	USH10M3H



Combiflange

Twist & Safe medium konvex transparent

Transparent, soft cover one side, tap valve.

Stoma opening	Plateau size	Article number
10 – 47 mm	55 mm	USH10M1FH
10 – 40 mm	45 mm	USH10M2FH
10 – 28 mm	35 mm	USH10M3FH



Combiflange

Stomocur® Safe konvex

7 mm
plateau height

Window, soft cover both sides, soft outlet, adhesive surface 116 mm (Ø), pouch size 155 x 210 mm, capacity 330 ml, package contents 10 pieces with 2 adapters.

Stoma opening	Plateau size	Article number
10 – 43 mm	50 mm	USH10V1F
15 mm		USH15V1F
20 mm		USH20V1F
25 mm		USH25V1F
30 mm		USH30V1F
35 mm		USH35V1F
10 – 35 mm	40 mm	USH10V2F
15 mm		USH15V2F
20 mm		USH20V2F
25 mm		USH25V2F
30 mm		USH30V2F
10 – 25 mm	30 mm	USH10V3F
15 mm		USH15V3F
20 mm		USH20V3F



Combiflange

Stomocur® Twist & Safe konvex

7 mm
plateau height

Window, soft cover both sides, tap valve, adhesive surface 116 mm (Ø), pouch size 155 x 210 mm, capacity 330 ml, package contents 10 pieces with 5 adapters & 5 drip caps.

Stoma opening	Plateau size	Article number
10 – 43 mm	50 mm	USH10V1H
15 mm		USH15V1H
20 mm		USH20V1H
25 mm		USH25V1H
30 mm		USH30V1H
35 mm		USH35V1H
10 – 35 mm	40 mm	USH10V2H
15 mm		USH15V2H
20 mm		USH20V2H
25 mm		USH25V2H
30 mm		USH30V2H
10 – 25 mm	30 mm	USH10V3H
15 mm		USH15V3H
20 mm		USH20V3H

Safe konvex transparent

Transparent, soft cover one side, soft outlet.

Stoma opening	Plateau size	Article number
10 – 43 mm	50 mm	USH10V1FT
10 – 35 mm	40 mm	USH10V2FT
10 – 25 mm	30 mm	USH10V3FT

Twist & Safe konvex transparent

Transparent, soft cover one side, tap valve.

Stoma opening	Plateau size	Article number
10 – 43 mm	50 mm	USH10V1HT
10 – 35 mm	40 mm	USH10V2HT
10 – 25 mm	30 mm	USH10V3HT



Konticur® Urostomy Bed Bag

Rotary tap drain for one-handed operation, kink-resistant tube, integrated adapter for connection to all urostomy pouches, package contents 5 pieces.

Tube length	Capacity	Article number
110 cm	2000 ml	KUBB2000



Stomocur® Uro Adapter

For connecting bed bags to the Stomocur® urostomy pouches, package contents 5 pieces.

Article number

USZ1



Stomocur® Drip Cap

For use with the Stomocur® 'Twist and Safe' tap, package contents 5 pieces.

Article number

USZ2





Stomocur®

Clic

Plate size 116 mm (Ø), package contents 5 pieces.



Locking ring	Stoma opening	Article number
45 mm	15 – 35 mm	CK4515
	20 mm	CK4520
	25 mm	CK4525
	30 mm	CK4530
57 mm	15 – 40 mm	CK5715
	35 mm	CK5735
	40 mm	CK5740
70 mm	15 – 50 mm	CK7015
	45 mm	CK7045
	50 mm	CK7050



Hydrocolloid



Stomocur®

Clic Alginat

Plate size 114 mm (Ø), package contents 5 pieces.



Locking ring	Stoma opening	Article number
45 mm	15 – 35 mm	CA4515
57 mm	15 – 40 mm	CA5715



Alginate



Stomocur®

Clic konvex

7 mm
plateau height

Plate size 116 mm (Ø), package contents 3 pieces.



Locking ring	Stoma opening	Plateau size	Article number
45 mm	15 – 35 mm	40 mm	CK4515V2
	20 mm		CK4520V2
	25 mm		CK4525V2
	30 mm		CK4530V2
	15 – 25 mm	30 mm	CK4515V3
57 mm	15 – 43 mm	50 mm	CK5715V1
	20 mm		CK5720V1
	25 mm		CK5725V1
	30 mm		CK5730V1
	35 mm		CK5735V1
	15 – 35 mm	40 mm	CK5715V2
	15 – 25 mm	30 mm	CK5715V3



Combiflange



Stomocur®

Clic Alginat konvex

7 mm
plateau height

Plate size 114 mm (Ø), package contents 3 pieces.



Locking ring	Stoma opening	Plateau size	Article number
45 mm	15 – 35 mm	40 mm	CA4515V2
	15 – 25 mm	30 mm	CA4515V3
57 mm	15 – 43 mm	50 mm	CA5715V1
	15 – 35 mm	40 mm	CA5715V2
	15 – 25 mm	30 mm	CA5715V3



Alginate



Stomocur®

Safe Clic

Window, soft cover both sides, soft outlet, pouch size 155 x 210 mm, capacity 330 ml, package contents 30 pieces with 5 adapters.

Locking ring	Article number
45 mm	CU45F
57 mm	CU57F

Safe Clic transparent

Transparent, soft cover one side, soft outlet.

Locking ring	Article number
45 mm	CU45FT
57 mm	CU57FT



Stomocur®

Twist & Safe Clic

Window, soft cover both sides, tap valve, pouch size 155 x 210 mm, capacity 330 ml, package contents 30 pieces with 5 adapters & 10 drip caps.

Locking ring	Article number
45 mm	CU45H
57 mm	CU57H

Twist & Safe Clic transparent

Transparent, soft cover one side, tap valve.

Locking ring	Article number
45 mm	CU45HT
57 mm	CU57HT



Special supply

Some situations call for solutions beyond everyday pouching systems. From neonatal care to complex drainage needs, specialised products help support safe and effective stoma management when standard options are not enough.

Our special supply range has been developed for medical professionals and offers specialised solutions that address these specific challenges, helping protect delicate skin, support clinical care and provide reassurance when it matters most.

Stomocur® Premature	68
Stomocur® Premature Combi	68
Stomocur® Drainage small - large	69
Stomocur® Drainage Extra "Access Port"	70

Paediatric care

- Two options for premature babies up to 900 g and infants up to 5000 g
- Gentle alginate infused hydrocolloid
- Soft non-woven cover on one side, transparent back for easy observation
- Closed pouch can be cut open and used as a drainable pouch with a clip

Drainage pouches

- Large skin protection plates combined with flexible, transparent pouches and integrated access port – ideal for easy observation
- High-output drainage ports with soft outlet
- A range of shapes and sizes to meet the patients' needs



Stomocur® Premature

For premature babies up to 900 g

Transparent, filter, soft cover one side, adhesive surface 40 mm, pouch size 50 x 120 mm, capacity 15 – 20 ml (with/without clip), package contents 20 pieces.

Stoma opening	Article number
3 mm	PA3003



Alginat



Stomocur® Premature Combi

For infants up to 5000 g

Transparent, soft cover one side, adhesive surface 40 mm, pouch size 73 x 220 mm, capacity 90 – 120 ml (with/without clip), package contents 30 pieces.

Stoma opening	Article number
3 mm	PA3203



Alginat

Stomocur® Drainage small - large

Transparent with printed scale, backstop, exhaust valve with adapter, package contents 30 pieces.



Hydrocolloid

Stoma opening	Capacity	Article number	Adhesive surface
3 mm	50 ml	DP2503XS	max. 54 mm
5 mm	150 ml	DP4005S	max. 92 mm
5 mm	300 ml	DP4005	max. 92 mm



Stomocur® Drainage Extra "Access Port"

Extra Small, opening 10 mm, transparent, adhesive surface 240 x 170 mm, pouch size 191 x 356 mm, package contents 1 piece.

Max. cutability	Access opening	Article number
90 x 180 mm	ø 120 mm	SDE2-XSA



Hydrocolloid



Stomocur® Drainage Extra "Access Port"

Small, opening 10 mm, transparent, adhesive surface 200 x 282 mm, pouch size 300 x 377 mm, package contents 1 piece.

Max. cutability	Access opening	Article number
240 x 155 mm	240 x 155 mm	SDE2-SA



Hydrocolloid



Stomocur® Drainage Extra "Access Port"

Medium, opening 10 mm, transparent, adhesive surface 305 x 134 mm, pouch size 177 x 430 mm, package contents 1 piece.

Max. cutability	Access opening	Article number
85 x 250 mm	85 x 250 mm	SDE2-MA



Hydrocolloid



Stomocur® Drainage Extra "Access Port"

Large, opening 10 mm, transparent, adhesive surface 262 x 221 mm, pouch size 250 x 442 mm, package contents 1 piece.

Max. cutability	Access opening	Article number
85 x 250 mm	160 x 200 mm	SDE2-LA



Hydrocolloid



Stomocur® Drainage Extra "Access Port"

With 2 drainage outlets

Extra Large, opening 10 mm, transparent, adhesive surface 450 x 300 mm, pouch size 350 x 500 mm, package contents 1 piece.

Max. cutability	Access opening	Article number
230 x 380 mm	ø 120 mm	SDE2-XXLA




Hydrocolloid

Accessories

Good stoma care extends beyond the pouching system itself. Small details – protecting the skin, improving adhesion or supporting comfortable wear – can make a meaningful difference in everyday life.

Our accessories bring together practical solutions designed to complement pouching systems and support peristomal skin health, helping create a reliable and comfortable foundation for daily stoma care.

Stomocur® Skin Protection Plates	74
Stomocur® Skin Protection Plates mouldable	74
Stomocur® Fixing Strips	74
Stomocur® Skin Protection Rings	75
Stomocur® Skin Protection Paste	75
Stomocur® Skin Protection Powder	75
Stomocur® Safety	76
Stomocur® EmplaSectal	76
Stomocur® Ostomy Belt	76
Stomocur® SGX	77
Stomocur® Scissors	77



Stomocur® Skin Protection Plates

With self-adhesive hydrocolloid skin protection compound, material thickness 0.6 mm, package contents 5 pieces.

Size	Article number
100 x 100 mm	HP101006
150 x 150 mm	HP151506
200 x 200 mm	HP202006
370 x 220 mm	HP372206



Stomocur® Skin Protection Plates mouldable

With self-adhesive hydrocolloid skin protection compound, can be individually cut and shaped, material thickness 2.4 mm, package contents 1 piece.

Size	Article number
220 x 220 mm	HP222224



Stomocur® Fixing Strips

Skin-friendly and highly elastic, extends the wearing time, prevents the baseplate from detaching and rolling up, more safety and comfort, package contents 20 pieces.

Article number
FS20



Stomocur® Skin Protection Rings

With hydrocolloid skin protection material to seal and to protect the skin around the stoma, available in two different outer diameters, material thickness 2.4 mm, package contents 30 pieces.

Outer Ø	Size	Article number
48	10 mm	48101
	15 mm	48151
	20 mm	48201
63	25 mm	63251
	30 mm	63301
	35 mm	63351
	40 mm	63401



Stomocur® Skin Protection Paste

For external use on the skin around the stoma, fills unevenness between the stoma and the skin protection plate, available with and without alcohol, package contents 1 piece.

Contents	Article number
56 g	CP2050 (with alcohol)
56 g	PASTEAO (without alcohol)



Stomocur® Skin Protection Powder

The powder absorbs the moisture of the skin around the stoma, ensuring better adhesion of the stoma supply, package contents 1 piece.

Contents	Article number
30 g	HP30



Stomocur® Safety

The skin care lotion is applied very thinly to the skin around the stoma to prevent skin irritation, package contents 1 piece.

Contents	Article number
100 ml	SSAFETY100



Stomocur® EmplaSectal

For easy and skin-friendly removal of the skin protection plates, handy spray, package contents 1 piece.

Contents	Article number
50 ml	EMPLASECTAL



Stomocur® Ostomy Belt

Ostomy belt for additional safety, available in two lengths, adjustable, width 3 cm, package contents 1 piece.

Lenght	Article number
100 cm	AB01
120 cm	AB02



Stomocur® SGX

Fast-acting, super-absorbent gelling agent, easy to use, helpful for ileostomists and colostomists with a high fluid output, reduces pouch swelling, pouch noise and odour, ideal for travel, package contents 30 pieces.

Article number
SGX30



Stomocur® Scissors

Scissors with rounded blades for cutting the baseplate openings, package contents 1 piece.

Article number
SSchere

Table of contents by article number

Article number	Opening	PZN	HIMI	Page
48101	10 mm	15865421	29.26.10.6011	75
48151	15 mm	15865438	29.26.10.6011	75
48201	20 mm	15865444	29.26.10.6011	75
63251	25 mm	15865450	29.26.10.6010	75
63301	30 mm	15865467	29.26.10.6010	75
63351	35 mm	15865473	29.26.10.6010	75
63401	40 mm	15865504	29.26.10.6010	75
AB01		6941287	29.26.11.0026	76
AB02		10525117	29.26.11.0027	76
CA4515V2	15 mm – 35 mm	18759027	29.26.06.0092	31, 47 & 63
CA4515V3	15 mm – 25 mm	18759033	29.26.06.0093	31, 47 & 63
CA5715V1	15 mm – 43 mm	18759056	29.26.06.0091	31, 47 & 63
CA5715V2	15 mm – 35 mm	18759062	29.26.06.0090	31, 47 & 63
CA5715V3	15 mm – 25 mm	18759079	29.26.06.0089	31, 47 & 63
CA4515	15 mm – 35 mm	18759085	29.26.05.0166	30, 46 & 62
CA5715	15 mm – 40 mm	18759091	29.26.05.0167	30, 46 & 62
CAA20V1B	20 mm – 43 mm	19536855	29.26.01.3085	28
CAA15V2B	15 mm – 35 mm	19536720	29.26.01.3086	28
CAA10V3B	10 mm – 25 mm	19536648	29.26.01.3087	28
CAA10M3B	10 mm – 28 mm	19536625	29.26.01.3089	27
CAA15M2B	15 mm – 40 mm	19536714	29.26.01.3088	27
CAA20M1B	20 mm – 47 mm	19536594	29.26.01.3090	27
CAH10V3B	10 mm – 25 mm	19536683	29.26.01.3081	27
CAH15V2B	15 mm – 35 mm	19536631	29.26.01.3080	27
CAH20V1B	20 mm – 43 mm	19536772	29.26.01.3079	27
CAH10M3B	10 mm – 28 mm	19536708	29.26.01.3084	26
CAH15M2B	15 mm – 40 mm	19536677	29.26.01.3083	26
CAH20M1B	20 mm – 47 mm	19536826	29.26.01.3082	26
CC45	45 mm	18238742	29.26.01.0029	32
CC57	57 mm	18238759	29.26.01.0037	32
CC70	70 mm	18238765	29.26.01.0038	32
CDP70T	70 mm	18884235	29.26.12.0018	45 & 49
CI45	45 mm	14367336	29.26.02.0105	48
CI45S	45 mm	16082980	29.26.02.0108	48

Article number	Opening	PZN	HIMI	Page
CI57	57 mm	14367365	29.26.02.0106	48
CI70	70 mm	16082974	29.26.02.0107	48
CK4515	15 mm – 35 mm	18238771	29.26.05.0156	30, 46 & 62
CK4515V2	15 mm – 35 mm	18238788	29.26.06.0062	31, 47 & 63
CK4515V3	15 mm – 25 mm	18238794	29.26.06.0061	31, 47 & 63
CK4520	20 mm	18238819	29.26.05.0156	30, 46 & 62
CK4520V2	20 mm	18238802	29.26.06.0062	31, 47 & 63
CK4525	25 mm	18239434	29.26.05.0156	30, 46 & 62
CK4525V2	25 mm	18238825	29.26.06.0062	31, 47 & 63
CK4530	30 mm	18238831	29.26.05.0156	30, 46 & 62
CK4530V2	30 mm	18238848	29.26.06.0062	31, 47 & 63
CK5715	15 mm – 40 mm	18238854	29.26.05.0157	30, 46 & 62
CK5715V1	15 mm – 43 mm	18238883	29.26.06.0063	31, 47 & 63
CK5715V2	15 mm – 35 mm	18238908	29.26.06.0064	31, 47 & 63
CK5715V3	15 mm – 25 mm	18238920	29.26.06.0086	31, 47 & 63
CK5720V1	20 mm	18238937	29.26.06.0063	31, 47 & 63
CK5725V1	25 mm	18238943	29.26.06.0063	31, 47 & 63
CK5730V1	30 mm	18238966	29.26.06.0063	31, 47 & 63
CK5735	35 mm	18238972	29.26.05.0157	30, 46 & 62
CK5735V1	35 mm	18238989	29.26.06.0063	31, 47 & 63
CK5740	40 mm	18238995	29.26.05.0157	30, 46 & 62
CK7015	15 mm – 50 mm	18239026	29.26.05.0158	30, 46 & 62
CK7045	45 mm	18239032	29.26.05.0158	30, 46 & 62
CK7050	50 mm	18239049	29.26.05.0158	30, 46 & 62
CP2050		7269751	29.26.10.0020	75
CSA10V3	10 mm – 25 mm	18172693	29.26.01.3051	29
CSA15V2	15 mm – 35 mm	18172894	29.26.01.3050	29
CSA15V3	15 mm	18192431	29.26.01.3051	29
CSA20	20 mm – 60 mm	18239055	29.26.01.2159	25
CSA20S	20 mm – 60 mm	14061583	29.26.01.2086	25
CSA20V1	20 mm – 43 mm	18198505	29.26.01.3052	29
CSA20V2	20 mm	18198511	29.26.01.3050	29
CSA20V3	20 mm	18198528	29.26.01.3051	29
CSA25	25 mm	18239061	29.26.01.2159	25
CSA25S	25 mm	14061614	29.26.01.2086	25
CSA25V1	25 mm	18198534	29.26.01.3052	29
CSA25V2	25 mm	18198540	29.26.01.3050	29
CSA30	30 mm	18239078	29.26.01.2159	25
CSA30S	30 mm	14061620	29.26.01.2086	25
CSA30V1	30 mm	18198824	29.26.01.3052	29
CSA30V2	30 mm	18198557	29.26.01.3050	29
CSA35	35 mm	18239084	29.26.01.2159	25

Article number	Opening	PZN	HIMI	Page
CSA35S	35 mm	14061643	29.26.01.2086	25
CSA35V1	35 mm	18198563	29.26.01.3052	29
CSA40	40 mm	18239090	29.26.01.2159	25
CSA40S	40 mm	14061666	29.26.01.2086	25
CSA40V1	40 mm	18198586	29.26.01.3052	29
CSC20	20 mm – 50 mm	14138185	29.26.04.0023	24
CSC20SNF	20 mm – 40 mm	18759116	29.26.04.0041	24
CSC25NF	25 mm – 50 mm	18759122	29.26.04.0042	24
CSC30LNF	30 mm – 60 mm	18759139	29.26.04.0043	24
CSG20	20 mm – 60 mm	17565545	29.26.01.2160	25
CSH10M3	10 mm – 28 mm	18198592	29.26.01.3044	26
CSH10V3	10 mm – 25 mm	18172902	29.26.01.3029	28
CSH10XS	10 mm – 50 mm	18172925	29.26.01.4013	24
CSH10XST	10 mm – 50 mm	16852857	29.26.01.4012	24
CSH15M2	15 mm – 40 mm	18198600	29.26.01.3043	26
CSH15V2	15 mm – 35 mm	18172919	29.26.01.3028	28
CSH15V3	15 mm	18192425	29.26.01.3029	28
CSH20	20 mm – 60 mm	18239109	29.26.01.2100	23
CSH20L	20 mm – 90 mm	18239150	29.26.01.2138	23
CSH20M1	20 mm – 47 mm	18198617	29.26.01.3040	26
CSH20M2	20 mm	18198623	29.26.01.3043	26
CSH20S	20 mm – 60 mm	12726268	29.26.01.2148	23
CSH20V0	20 mm – 57 mm	18198646	29.26.01.3049	28
CSH20V1	20 mm – 43 mm	18198675	29.26.01.3027	28
CSH20V2	20 mm	18198681	29.26.01.3028	28
CSH20V3	20 mm	18198698	29.26.01.3029	28
CSH25	25 mm	18239115	29.26.01.2100	23
CSH25L	25 mm	18239167	29.26.01.2138	23
CSH25M1	25 mm	18198712	29.26.01.3040	26
CSH25M2	25 mm	18198729	29.26.01.3043	26
CSH25S	25 mm	12726251	29.26.01.2148	23
CSH25V1	25 mm	18198735	29.26.01.3027	28
CSH25V2	25 mm	18198741	29.26.01.3028	28
CSH30	30 mm	18239121	29.26.01.2100	23
CSH30L	30 mm	18239173	29.26.01.2138	23
CSH30M1	30 mm	18198758	29.26.01.3040	26
CSH30S	30 mm	12726239	29.26.01.2148	23
CSH30V1	30 mm	18198764	29.26.01.3027	28
CSH30V2	30 mm	18198770	29.26.01.3028	28
CSH35	35 mm	18239138	29.26.01.2100	23
CSH35L	35 mm	18239196	29.26.01.2138	22
CSH35M1	35 mm	18198787	29.26.01.3040	26

Article number	Opening	PZN	HIMI	Page
CSH35S	35 mm	12726222	29.26.01.2148	23
CSH35V1	35 mm	18198793	29.26.01.3027	28
CSH40	40 mm	18239144	29.26.01.2100	23
CSH40L	40 mm	18239204	29.26.01.2138	22
CSH40M1	40 mm	18198801	29.26.01.3040	26
CSH40S	40 mm	12726245	29.26.01.2148	23
CSH40V1	40 mm	18198818	29.26.01.3027	28
CSIL15B	15 – 60 mm	20546497	29.26.01.2900	22
CU45F	45 mm	17206869	29.26.03.0060	64
CU45FT	45 mm	17206823	29.26.03.0061	64
CU45H	45 mm	17554866	29.26.03.0064	64
CU45HT	45 mm	17554889	29.26.03.0063	64
CU57F	57 mm	17206875	29.26.03.0059	64
CU57FT	57 mm	17206852	29.26.03.0062	64
CU57H	57 mm	17554872	29.26.03.0066	64
CU57HT	57 mm	17554895	29.26.03.0065	64
DP10M3	10 mm – 28 mm	18884287	29.26.12.3009	42 & 45
DP10V3	10 mm – 25 mm	18884318	29.26.12.3003	45
DP15M2	15 mm – 40 mm	18884270	29.26.12.3010	42 & 45
DP15V2	15 mm – 35 mm	18884301	29.26.12.3002	45
DP20M1	20 mm – 47 mm	18884264	29.26.12.3008	42 & 45
DP20V1	20 mm – 43 mm	18884293	29.26.12.3001	45
DP2503XS	3 mm	18239210	29.26.12.2015	69
DP4005	5 mm	18239227	29.26.12.2012	69
DP4005S	5 mm	18239233	29.26.12.2014	69
DP5020	20 mm – 60 mm	18884324	29.26.12.2013	40 & 45
DPB2000		18827467	29.26.11.6900	45
EMPLASECTAL		200420	29.26.10.4020	76
FS20		16622583	29.26.11.6173	74
HP101006	100 mm x 100 mm	10627941	29.26.07.0045	74
HP151506	150 mm x 150 mm	18239279	29.26.07.1024	74
HP202006	200 mm x 200 mm	10627935	29.26.07.2016	74
HP222224	220 mm x 220 mm	18239316	29.26.07.2015	74
HP30		4572710	29.26.10.1009	75
HP372206	370 mm x 220 mm	4877390	29.26.07.2014	74
IAA10M3B	10 mm – 28 mm	19536766	29.26.02.3125	42
IAA15M2B	15 mm – 40 mm	19536654	29.26.02.3124	42
IAA15M1B	15 mm – 47 mm	19536789	29.26.02.3123	42
IAA10V3B	10 mm – 25 mm	19536795	29.26.02.3120	44
IAA15V1B	15 mm – 43 mm	19536832	29.26.02.3118	44
IAA15V2B	15 mm – 35 mm	19536737	29.26.02.3119	44
IAH10M3B	10 mm – 28 mm	19536602	29.26.02.3126	41

Article number	Opening	PZN	HIMI	Page
IAH15M1B	15 mm – 47 mm	19536660	29.26.02.3121	41
IAH15M2B	15 mm – 40 mm	19536849	29.26.02.3122	41
IAH15V1B	15 mm – 43 mm	19536861	29.26.02.3115	43
IAH15V2B	15 mm – 35 mm	19536619	29.26.02.3116	43
IAH10V3B	10 mm – 25 mm	19536803	29.26.02.3117	43
IPA10V3	10 mm – 25 mm	18192514	29.26.02.3081	44
IPA15	15 mm – 60 mm	14419990	29.26.02.2172	39
IPA15V1	15 mm – 43 mm	18192744	29.26.02.3079	44
IPA15V2	15 mm – 35 mm	18192508	29.26.02.3080	44
IPA20	20 mm	14363491	29.26.02.2172	39
IPA20V1	20 mm	14367129	29.26.02.3079	44
IPA25	25 mm	14363516	29.26.02.2172	39
IPA25V1	25 mm	14367135	29.26.02.3079	44
IPA30	30 mm	14363522	29.26.02.2172	39
IPA30V1	30 mm	14367141	29.26.02.3079	44
IPA35	35 mm	14363539	29.26.02.2172	39
IPA35V1	35 mm	16128794	29.26.02.3079	44
IPA40	40 mm	14363640	29.26.02.2172	39
IPA40V1	40 mm	18192477	29.26.02.3079	44
IPG20	20 mm – 60 mm	17565551	29.26.02.2181	40
IPH10M3	10 mm – 28 mm	14360274	29.26.02.3086	41
IPH10M3T	10 mm – 28 mm	17551431	29.26.02.3093	41
IPH10V3	10 mm – 25 mm	18192566	29.26.02.3088	43
IPH15	15 mm – 60 mm	14419978	29.26.02.2179	38
IPH15M1	15 mm – 47 mm	18192454	29.26.02.3084	41
IPH15M1T	15 mm – 47 mm	18192460	29.26.02.3092	41
IPH15M2	15 mm – 40 mm	14360245	29.26.02.3085	41
IPH15M2T	15 mm – 40 mm	17551419	29.26.02.3091	41
IPH15T	15 mm – 60 mm	16134174	29.26.02.2173	38
IPH15S	15 mm – 60 mm	17551336	29.26.02.2186	39
IPH15V1	15 mm – 43 mm	14334638	29.26.02.3083	43
IPH15V2	15 mm – 35 mm	18192543	29.26.02.3087	43
IPH15V3	15 mm	18192448	29.26.02.3088	43
IPH20	20 mm	14330362	29.26.02.2179	38
IPH20M1	20 mm	14360239	29.26.02.3084	41
IPH20M2	20 mm	17617408	29.26.02.3085	41
IPH20V0	20 mm – 57 mm	14334710	29.26.02.3082	43
IPH20V1	20 mm	14334644	29.26.02.3083	43
IPH20V2	20 mm	14360185	29.26.02.3087	43
IPH20V3	20 mm	14360222	29.26.02.3088	43
IPH2070	20 mm – 70 mm	18825712	29.26.02.2187	38
IPH25	25 mm	14330379	29.26.02.2179	38

Article number	Opening	PZN	HIMI	Page
IPH25M1	25 mm	17617377	29.26.02.3084	41
IPH25M2	25 mm	17558410	29.26.02.3085	41
IPH25V1	25 mm	14334650	29.26.02.3083	43
IPH25V2	25 mm	14360191	29.26.02.3087	43
IPH30	30 mm	14330385	29.26.02.2179	38
IPH30M1	30 mm	17558350	29.26.02.3084	41
IPH30V1	30 mm	14334667	29.26.02.3083	43
IPH30V2	30 mm	14360216	29.26.02.3087	43
IPH35	35 mm	14330416	29.26.02.2179	38
IPH35M1	35 mm	17558373	29.26.02.3084	41
IPH35V1	35 mm	14334673	29.26.02.3083	43
IPH40	40 mm	14330422	29.26.02.2179	38
IPH40M1	40 mm	17558404	29.26.02.3084	41
IPH40V1	40 mm	14334704	29.26.02.3083	43
ISIL15B	15 – 60 mm	20546505	29.26.02.2900	38
KUBB2000		18095774	29.26.11.6900	60
PA3003	3 mm	20322869		68
PA3203	3 mm	20322875		68
PASTEOA		14055921	29.26.10.0038	75
SDE2–LA	160 mm x 200 mm	11549108	29.26.12.2009	71
SDE2–MA	85 mm x 250 mm	11549143	29.26.12.2007	71
SDE2–SA	240 mm x 155 mm	11549172	29.26.12.2006	70
SDE2–XSA	90 mm x 180 mm	11549203	29.26.12.2011	70
SDE2–XXLA	230 mm x 380 mm	11549120	29.26.12.2010	71
SGX30		17511897		77
SSAFETY100		4751772	29.26.10.3019	76
SSchere				77
USG15FT	15 mm – 60 mm	17255980	29.26.03.2061	56
USH10F	10 mm – 60 mm	18192572	29.26.03.2053	54
USH10FT	10 mm – 60 mm	18192589	29.26.03.2060	54
USH10H	15 mm – 60 mm	18192595	29.26.03.2065	54
USH10HT	10 mm – 60 mm	18192603	29.26.03.2066	54
USH10M1F	10 mm – 47 mm	18192626	29.26.03.3034	56
USH10M1FT	10 mm – 47 mm	18192655	29.26.03.3044	56
USH10M1H	10 mm – 47 mm	18192684	29.26.03.3061	57
USH10M1HT	10 mm – 47 mm	18192715	29.26.03.3060	57
USH10M2F	10 mm – 40 mm	18192649	29.26.03.3033	56
USH10M2FT	10 mm – 40 mm	18192678	29.26.03.3049	56
USH10M2H	10 mm – 40 mm	18192709	29.26.03.3063	57
USH10M2HT	10 mm – 40 mm	18192738	29.26.03.3052	57
USH10M3F	10 mm – 28 mm	13867002	29.26.03.3032	56
USH10M3FT	10 mm – 28 mm	16622459	29.26.03.3046	56

Article number	Opening	PZN	HIMI	Page
USH10M3H	10 mm – 28 mm	17557988	29.26.03.3062	57
USH10M3HT	10 mm – 28 mm	17557965	29.26.03.3056	57
USH10V1H	10 mm – 43 mm	17552867	29.26.03.3053	59
USH10V1HT	10 mm – 43 mm	17552927	29.26.03.3059	59
USH10V1F	10 mm – 43 mm	13866793	29.26.03.3035	58
USH10V1FT	10 mm – 43 mm	16622465	29.26.03.3045	58
USH10V2F	10 mm – 35 mm	13866801	29.26.03.3036	58
USH10V2FT	10 mm – 35 mm	16622376	29.26.03.3048	58
USH10V2H	10 mm – 35 mm	17552904	29.26.03.3055	59
USH10V2HT	10 mm – 35 mm	17552962	29.26.03.3058	59
USH10V3F	10 mm – 25 mm	13866818	29.26.03.3037	58
USH10V3FT	10 mm – 25 mm	16622554	29.26.03.3047	58
USH10V3H	10 mm – 25 mm	17552991	29.26.03.3054	59
USH10V3HT	10 mm – 25 mm	17553045	29.26.03.3057	59
USH15F	15 mm	17160535	29.26.03.2053	54
USH15H	15 mm	17552809	29.26.03.2065	54
USH15M1F	15 mm	18192632	29.26.03.3034	56
USH15M1H	15 mm	18192690	29.26.03.3061	57
USH15M2F	15 mm	13866994	29.26.03.3033	56
USH15M2H	15 mm	17557942	29.26.03.3063	57
USH15V1F	15 mm	13866824	29.26.03.3035	58
USH15V1H	15 mm	17757084	29.26.03.3053	59
USH15V2F	15 mm	13866830	29.26.03.3036	58
USH15V2H	15 mm	17757109	29.26.03.3055	59
USH15V3F	15 mm	13866847	29.26.03.3037	58
USH15V3H	15 mm	17553016	29.26.03.3054	59
USH20F	20 mm	17160541	29.26.03.2053	54
USH20H	20 mm	17552815	29.26.03.2065	54
USH20M1F	20 mm	13866988	29.26.03.3034	56
USH20M1H	20 mm	17557899	29.26.03.3061	57
USH20M2F	20 mm	17622355	29.26.03.3033	56
USH20M2H	20 mm	17757196	29.26.03.3063	57
USH20V1F	20 mm	13866876	29.26.03.3035	58
USH20V1H	20 mm	17757078	29.26.03.3053	59
USH20V2F	20 mm	13866882	29.26.03.3036	58
USH20V3F	20 mm	13866899	29.26.03.3037	58
USH20V2H	20 mm	17757115	29.26.03.3055	59
USH20V3H	20 mm	17553022	29.26.03.3054	59
USH25F	25 mm	17160558	29.26.03.2053	54
USH25H	25 mm	17552821	29.26.03.2065	54
USH25M1F	25 mm	17617466	29.26.03.3034	56
USH25M1H	25 mm	17757173	29.26.03.3061	57

Article number	Opening	PZN	HIMI	Page
USH25M2F	25 mm	17557824	29.26.03.3033	56
USH25M2H	25 mm	17557971	29.26.03.3063	57
USH25V1F	25 mm	13866913	29.26.03.3035	58
USH25V1H	25 mm	17757090	29.26.03.3053	59
USH25V2F	25 mm	13866936	29.26.03.3036	58
USH25V2H	25 mm	17552910	29.26.03.3055	59
USH30F	30 mm	17160564	29.26.03.2053	54
USH30H	30 mm	17552838	29.26.03.2065	54
USH30M1F	30 mm	17557770	29.26.03.3034	56
USH30M1H	30 mm	17557907	29.26.03.3061	57
USH30V1F	30 mm	13866942	29.26.03.3035	58
USH30V1H	30 mm	17552873	29.26.03.3053	59
USH30V2F	30 mm	13866959	29.26.03.3036	58
USH30V2H	30 mm	17552985	29.26.03.3055	59
USH35M1F	35 mm	17557787	29.26.03.3034	56
USH35M1H	35 mm	17557913	29.26.03.3061	57
USH35V1F	35 mm	13866965	29.26.03.3035	58
USH35V1H	35 mm	17552896	29.26.03.3053	59
USM13F	13 mm – 50 mm	17160386	29.26.03.1014	55
USM13FT	13 mm – 50 mm	17206504	29.26.03.1015	55
USM13H	13 mm – 50 mm	17846770	29.26.03.1018	55
USM13HT	13 mm – 50 mm	17846787	29.26.03.1017	55
USZ1		13867019	29.26.10.9000	60
USZ2		17931464	29.26.10.9000	60



Here for what matters

Unora GmbH
Am Studio 16
12489 Berlin
Germany

mail@unora.de
www.unora.de

© Unora GmbH 2026

unora