

Product guide

unoro

2026

Here for what matters



Welcome to Unora

We're here to redefine what's possible in stoma care – for people living with change and for you, who share our commitment to helping them live the lives they deserve.

Our holistic and flexible portfolio of products and accessories is designed with protection and comfort in mind.

Combining innovative materials and decades of experience, our five material options make it easy to shape care to each person – whatever their needs and wherever they are on their journey.

And because care goes beyond the product itself, we're here for the quiet victories, the hard moments and everything in between. From the initial conversations to life after surgery, we're the partner and lifeline, ready with guidance and support that helps everyone feel confident every step of the way.

With Unora, you can focus on what matters most – delivering the best care to more people.

We're here for what matters.

New product range

Explore the future of stoma care

Introducing the first ostomy bag with skin-friendly silicone

Silkura® is the only silicone based stoma bag range and is designed to offer breathable adhesion, comfort and secure wear. Its silicone barrier allows airflow to help support the skin's natural moisture balance, while adapting closely to the skin to help create a secure seal during everyday activity. The flexible flange moves naturally with the body for a light, barely noticeable feel.



Closed and drainable bags in a discreet black design with soft textile cover. Flat baseplate, suitable for ileostomy and colostomy.

Find our Silkura® pouches on pages 22 and 34.



Alginate



Alginate infused hydrocolloid designed to help keep skin drier and more comfortable

Unora's Alginate baseplate is formulated with 15% alginate, a naturally derived ingredient known for its moisture absorbing, soothing properties. As it absorbs fluid, it forms a soft, elastic gel that adapts to the skin surface – supporting comfort, adhesion and peristomal skin integrity.

With up to 1.5 times greater moisture absorption than standard hydrocolloids*, the Alginate range is ideal for:

- Moist or perspiration prone peristomal skin
- Skin prone to allergies
- Ulcerated or sore peristomal skin

Available options

Closed and drainable bags in beige and selected black styles. Flat, low profile (3.5 mm) and soft convex (7 mm) baseplates. Suitable for ileostomy and colostomy.

* Based on standardised laboratory testing showing up to 1.5 times higher fluid absorption compared to our standard hydrocolloid materials.

Here for what matters



Unora Home Delivery provides patients with a reliable and compassionate home delivery service for all stoma prescription products. As one of the UK's leading services, we supply the full Unora range alongside products from all major manufacturers, ensuring choice, continuity and confidence.

Unora Home Delivery®

Supporting patients

Our service has been designed around what matters most to people living with a stoma – ease, dignity and dependable support.

We offer:

- **Free, discreet home delivery** anywhere in the UK
- **Electronic Prescription Service (EPS)** or GP collection to streamline ordering
- Access to **Unora Clinical Nurse Specialists**, highly experienced in stoma care
- **Complimentary items** such as wipes, disposal bags and practical comfort essentials
- A dedicated **Customer Care Team** offering impartial support Mon–Fri

Virtual nursing service

For patients needing additional guidance, our Virtual Nursing Service offers secure video or telephone consultations with a Unora stoma nurse – ideal for follow up support, product advice, trouble-shooting and wellbeing conversations. It complements, not replaces, local NHS care.

Why Unora?

Built on decades of clinical expertise and a bold commitment to making care more human, Unora stands for clarity, dignity and genuine presence in a system that can often feel impersonal. We're here for what matters – to you, and to your patients.

Simple registration

Healthcare professionals can refer patients directly, or patients can register themselves by phone or email.

They will need:

- Contact details
- GP information
- Medical exemption certificate number (if applicable)

Freephone 0800 592 786
Email priority@unorahomedelivery.co.uk

Contents

Introduction

Skin protection materials	12
Convexities	14
Product features	16

Products

Colostomy	18
Ileostomy	30
Urostomy	42
Wound management	54
Accessories	60



Skin protection materials

Protecting peristomal skin is central to effective stoma care. Our portfolio includes a range of skin protection materials designed to support different skin needs and wearing preferences, making it easier to share care around the individual.



Silicone

Breathable silicone skin protection that combines gentle adhesion with a secure seal.

- Breathable barrier that helps maintain the skin's natural moisture balance
- Adapts closely to the skin's surface to help create a secure seal
- Moves naturally with the body for a light, barely noticeable feel



Alginate

Alginate-infused skin protection designed to absorb up to 1.5 times more fluid than standard hydrocolloids* for improved moisture management.

- Contains 15% alginate derived from natural seaweed
- Forms an elastic gel that adapts to any skin surface
- Helps maintain a balanced peristomal environment and supports comfort and extended wear time

* Based on standardised laboratory testing showing up to 1.5 times higher absorption compared with standard hydrocolloid materials. Data on file.



Primo hydrocolloid

Dual-layer hydrocolloid technology engineered to enhance adhesion and security while remaining flexible and comfortable.

- Two hydrocolloid layers bonded together to increase absorptive capacity
- 0.4 mm thicker structure compared to standard hydrocolloid for enhanced adhesion
- Thin, tapered border designed to adapt to uneven skin contours



Microskin

Ultra-thin, flexible skin protection that adapts closely to the body and supports natural movement.

- Ultra-thin foil with a thickness of just 0.02 mm thick for improved discretion
- Highly flexible design that adapts to body contours and movement
- Suitable for active lifestyles



Hydrocolloid


Used in our wound drainage products, this hydrocolloid provides skin protection with reliable adhesion and a secure seal to fit everyday needs.

- Established hydrocolloid formulation for dependable adhesion
- Absorbs moisture around the stoma to help maintain a secure seal
- Proven material technology widely used in everyday ostomy care

Convexities


Every stoma and body contour is different. When additional support is needed to help the stoma protrude and maintain a secure seal, convex baseplates can help create a more reliable fit.

Our convex portfolio offers different convexity heights and plateau sizes, making it easier to match the right level of support to each person's needs.



Low Profile Soft Convex

Low Profile Soft Convex baseplates provide gentle support while maintaining flexibility and comfort. With a convex height of 3.5 mm, they help create a more secure fit while adapting naturally to the body's contours. Low Profile Soft Convex baseplates are available in three different plateau sizes.



Soft Convex

Soft Convex baseplates provide a deeper level of support when additional assistance is needed to help the stoma protrude and maintain a secure seal. With a convex height of 7 mm, they are designed to support more challenging stoma profiles while maintaining comfort and flexibility. Soft Convex baseplates are available in up to four different plateau sizes.

Available in Primo Hydrocolloid and Alginate materials, our Low Profile Soft Convex and Soft Convex options combine reliable adhesion with soft, flexible construction designed for comfortable everyday wear.


Product features

Our pouches are developed and manufactured in Berlin using modern production technologies and high-quality materials.


They are characterised by innovative and unique properties that meet the medical requirements of healthcare professionals and users.

Pouch shapes


Symmetric



Asymmetric




Outlets on high output pouches




Cover material

Non-woven cover



Textile cover




Filter

- Colostomy – Activated carbon filter
- Ileostomy – MPC filter
MPC filter technology (Moisture Protect Composition)
= additional protective layer that prevents liquid absorption.

Window

All of our bags feature a window for easy positioning and observation. Our clear bags without cover/ window feature a crystal clear film, again for easy observation.



Colostomy

Unora colostomy solutions combine reliable filters, comfortable materials and skin-friendly barriers to help protect the peristomal area and maintain confidence in daily wear.

The Unora colostomy portfolio spans four skin protection materials that make it easier to shape care around individual anatomy, skin condition and lifestyle – whether supporting established routines or helping someone regain confidence after surgery.



Silicone

Breathable silicone skin protection that combines gentle adhesion with a secure seal.



Alginate

Alginate-infused skin protection designed to absorb up to 1.5 times more fluid than standard hydrocolloids* for improved moisture management.



Primo hydrocolloid

Dual-layer hydrocolloid technology engineered to enhance adhesion and security while remaining flexible and comfortable.



Microskin

Ultra-thin, flexible skin protection that adapts closely to the body and supports natural movement.

* Based on standardised laboratory testing showing up to 1.5 times higher absorption compared with standard hydrocolloid materials. Data on file.

Product overview



Silicone

Silkura® 22



Alginate

Closed Standard Alginate 22

Closed Mini Alginate 23

Closed Maxi Alginate 23

Closed Standard Alginate
Low Profile Soft Convex 24

Closed Standard Alginate
Soft Convex 24



Primo hydrocolloid

Closed Mini Primo Hydrocolloid 25

Closed Standard Primo Hydrocolloid 25

Closed Maxi Primo Hydrocolloid 26

Closed Standard Primo Hydrocolloid
Low Profile Soft Convex 26

Closed Standard Primo Hydrocolloid
Soft Convex 27



Microskin

Closed Standard Microskin 28

Closed Plus Microskin 28



Silkura®

New for 2026

Window, filter, ultra-soft, opaque and water-repellent cover on both sides, adhesive surface max. 130 mm, pouch size 154 x 202 mm, capacity 330 ml.
Pack size: 30



Silicone

Hole size	Window
Cut to fit 15 – 60 mm	TR322C



Closed Standard Alginate

AHP5C20 - AHP5C60 Window, opaque, filter, soft cover both sides.
Pack Size: 30



Alginate

Hole size	Window
Cut to fit 20 to 60 mm	AHP5C20
Pre-cut to 25 mm	AHP5C25
Pre-cut to 30 mm	AHP5C30
Pre-cut to 35 mm	AHP5C35
Pre-cut to 40 mm	AHP5C40
Pre-cut to 45 mm	AHP5C45
Pre-cut to 50 mm	AHP5C50
Pre-cut to 55 mm	AHP5C55
Pre-cut to 60 mm	AHP5C60



Closed Mini Alginate

AHP11C10 Window, opaque, filter, soft cover both sides.
Pack Size: 30



Alginate

Hole size	Window
Cut to fit 10 to 50 mm	AHP11C10



Closed Maxi Alginate

AHP8C20 Window, opaque, filter, soft cover both sides.
AHP7C20C Clear, soft cover one side.
Pack Size: 30



Alginate

Hole size	Window
Cut to fit 20 to 90 mm	AHP8C20



Closed Standard Alginate Low Profile Soft Convex

3.5 mm plateau height

LPSCA28-0410KBK - LPSCA47-0445K Window, opaque, filter, soft cover both sides.
Pack Size: 20

Hole size	Plateau	Window
● Cut to fit 10 to 28 mm	28 mm	LPSCA28-0410KBK
● Cut to fit 20 to 28 mm	28 mm	LPSCA28-0420K
● Pre-cut to 25 mm	28 mm	LPSCA28-0425K
● Cut to fit 15 to 40 mm	40 mm	LPSCA40-0415KBK
● Cut to fit 20 to 40 mm	40 mm	LPSCA40-0420K
● Pre-cut to 30 mm	40 mm	LPSCA40-0430K
● Pre-cut to 35 mm	40 mm	LPSCA40-0435K
● Cut to fit 20 to 47 mm	47 mm	LPSCA47-0420KBK
● Cut to fit 20 to 47 mm	47 mm	LPSCA47-0420K
● Pre-cut to 40 mm	47 mm	LPSCA47-0440K
● Pre-cut to 45 mm	47 mm	LPSCA47-0445K



Alginate

● Black cover due Autumn 2026 ● Beige cover



Closed Standard Alginate Soft Convex

7 mm plateau height

SCA25-0410KBK - SCA45-0440K Window, opaque filter, soft cover both sides.
Pack Size: 20

Hole size	Plateau	Window
● Cut to fit 10 to 24 mm	25 mm	SCA25-0410KBK
● Cut to fit 20 to 24 mm	25 mm	SCA25-0420K
● Cut to fit 15 to 34 mm	35 mm	SCA35-0415KBK
● Cut to fit 20 to 34 mm	35 mm	SCA35-0420K
● Pre-cut to 25 mm	35 mm	SCA35-0425K
● Pre-cut to 30 mm	35 mm	SCA35-0430K
● Cut to fit 20 to 44 mm	45 mm	SCA45-0420KBK
● Cut to fit 20 to 44 mm	45 mm	SCA45-0420K
● Pre-cut to 35 mm	45 mm	SCA45-0435K
● Pre-cut to 40 mm	45 mm	SCA45-0440K



Alginate

● Black cover due Autumn 2026 ● Beige cover



Closed Mini Primo Hydrocolloid

HP11C10 Window, opaque, filter, soft cover both sides.
Pack Size: 30

Hole size	Window
Cut to fit 10 to 50 mm	HP11C10



Primo hydrocolloid




Closed Standard Primo Hydrocolloid

HP5C20 - HP5C60 Window, opaque, filter, soft cover both sides.
HP5C20C Clear, filter, soft cover one side.
Pack Size: 30

Hole size	Window	Clear
Cut to fit 20 to 60 mm	HP5C20	HP5C20C
Pre-cut to 25 mm	HP5C25	
Pre-cut to 30 mm	HP5C30	
Pre-cut to 35 mm	HP5C35	
Pre-cut to 40 mm	HP5C40	
Pre-cut to 45 mm	HP5C45	
Pre-cut to 50 mm	HP5C50	
Pre-cut to 55 mm	HP5C55	
Pre-cut to 60 mm	HP5C60	



Primo hydrocolloid



Closed Maxi Primo Hydrocolloid

HP8C20 Window, opaque, filter, soft cover both sides.


HP7C20C Clear, filter, soft cover one side.

Pack Size: 30

Hole size	Window	Clear
Cut to fit 20 to 90 mm	HP8C20	HP7C20C



Primo hydrocolloid



Closed Standard Primo Hydrocolloid Low Profile Soft Convex

3.5 mm
plateau height

LPSC28-0410KBK - LPSC47-0445K Window,

opaque, filter, soft cover both sides.


Pack Size: 20

Hole size	Plateau	Window
● Cut to fit 10 to 28 mm	28 mm	LPSC28-0410KBK
● Cut to fit 20 to 28 mm	28 mm	LPSC28-0420K
● Pre-cut to 25 mm	28 mm	LPSC28-0425K
● Cut to fit 15 to 40 mm	40 mm	LPSC40-0415KBK
● Cut to fit 20 to 40 mm	40 mm	LPSC40-0420K
● Pre-cut to 30 mm	40 mm	LPSC40-0430K
● Pre-cut to 35 mm	40 mm	LPSC40-0435K
● Cut to fit 20 to 47 mm	47 mm	LPSC47-0420KBK
● Cut to fit 20 to 47 mm	47 mm	LPSC47-0420K
● Pre-cut to 40 mm	47 mm	LPSC47-0440K
● Pre-cut to 45 mm	47 mm	LPSC47-0445K

● Black cover due Autumn 2026 ● Beige cover



Primo hydrocolloid



Closed Standard Primo Hydrocolloid Soft Convex

7 mm
plateau height

SC25-0410KBK - SC57-0420K Window,

opaque filter, soft cover both sides.

Pack Size: 20

Hole size	Plateau	Window
● Cut to fit 10 to 24 mm	25 mm	SC25-0410KBK
● Cut to fit 20 to 24 mm	25 mm	SC25-0420K
● Cut to fit 15 to 34 mm	35 mm	SC35-0415KBK
● Cut to fit 20 to 34 mm	35 mm	SC35-0420K
● Pre-cut to 25 mm	35 mm	SC35-0425K
● Pre-cut to 30 mm	35 mm	SC35-0430K
● Cut to fit 20 to 44 mm	45 mm	SC45-0420KBK
● Cut to fit 20 to 44 mm	45 mm	SC45-0420K
● Pre-cut to 35 mm	45 mm	SC45-0435K
● Pre-cut to 40 mm	45 mm	SC45-0440K



Primo hydrocolloid

Pack Size: 10

Hole size	Plateau	Window
● Cut to fit 20 to 57 mm	57 mm	SC57-0420K

● Black cover due Autumn 2026 ● Beige cover



Closed Standard Microskin

M-0320K - M-0340K Window, filter,
soft cover both sides.

Pack Size: 30

Hole size	Window
Cut to fit 20 to 50 mm	M-0320K
Pre-cut to 25 mm	M-0325K
Pre-cut to 30 mm	M-0330K
Pre-cut to 35 mm	M-0335K
Pre-cut to 40 mm	M-0340K



Microskin



Closed Plus Microskin

M-0420K - M-0440K Window, filter,
soft cover both sides.

Pack Size: 30

Hole size	Window
Cut to fit 20 to 50 mm	M-0420K
Pre-cut to 25 mm	M-0425K
Pre-cut to 30 mm	M-0430K
Pre-cut to 35 mm	M-0435K
Pre-cut to 40 mm	M-0440K



Microskin



Ileostomy

Ileostomy care often places greater demands on both products and skin. Higher output, moisture and frequent emptying mean that reliable sealing, skin protection and ease of use are essential.

Our ileostomy solutions are designed to support throughout the day. Advanced filter technology, secure closures and skin-friendly materials work together to protect the peristomal skin and help reduce leakage risk.

Our ileostomy portfolio offers a range of skin protection materials and pouch options to help adapt care to different anatomies, skin needs and daily routines.



Silicone

Breathable silicone skin protection that combines gentle adhesion with a secure seal.



Alginate

Alginate-infused skin protection designed to absorb up to 1.5 times more fluid than standard hydrocolloids* for improved moisture management.



Primo hydrocolloid

Dual-layer hydrocolloid technology engineered to enhance adhesion and security while remaining flexible and comfortable.



Microskin

Ultra-thin, flexible skin protection that adapts closely to the body and supports natural movement.

* Based on standardised laboratory testing showing up to 1.5 times higher absorption compared with standard hydrocolloid materials. Data on file.

Product overview



Silicone

Silkura® 34



Alginate

Drainable Standard Alginate 34

Drainable Midi Alginate 35

Drainable Standard Alginate
Low Profile Soft Convex 35

Drainable Standard Alginate Soft Convex 36



Primo hydrocolloid

Drainable Midi Primo Hydrocolloid 36

Drainable Standard Primo Hydrocolloid 37

Drainable Standard Primo Hydrocolloid
High Output 37

Drainable Standard Primo Hydrocolloid
High Output Soft Convex 38

Drainable Midi Primo Hydrocolloid
Low Profile Soft Convex 38

Drainable Standard Primo Hydrocolloid
Low Profile Soft Convex 39

Drainable Midi Primo Hydrocolloid
Soft Convex 40

Drainable Standard Primo Hydrocolloid
Soft Convex 41



Microskin

Drainable Standard Microskin 41



Silkura®

New for 2026

Window, filter, ultra-soft, opaque and water-repellent cover on both sides, adhesive surface max. 130 mm, pouch size 155 x 215 mm (fully closed), capacity 460 ml. Pack size: 30

Hole size	Window
Cut to fit 15 – 60 mm	TR325C



Silicone



Drainable Standard Alginate

AV-4120K - AV-4160K Window, filter, soft cover both sides, soft touch closure. Pack Size: 30

Hole size	Window
Cut to fit 20 to 60 mm	AV-4120K
Pre-cut to 25 mm	AV-4125K
Pre-cut to 30 mm	AV-4130K
Pre-cut to 35 mm	AV-4135K
Pre-cut to 40 mm	AV-4140K
Pre-cut to 45 mm	AV-4145K
Pre-cut to 50 mm	AV-4150K
Pre-cut to 55 mm	AV-4155K
Pre-cut to 60 mm	AV-4160K



Alginate



Drainable Midi Alginate

AV-4320K Window, filter, soft cover both sides, soft touch closure. Pack Size: 30

Hole size	Window Midi
Cut to fit 20 to 60 mm	AV-4320K



Alginate



Drainable Standard Alginate Low Profile Soft Convex

3.5 mm plateau height

LPSCA28-4110KVBK - LPSCA47-4145KV Window, opaque, filter, soft cover both sides, soft touch closure. Pack Size: 20

Hole size	Plateau	Window
● Cut to fit 10 to 28 mm	28 mm	LPSCA28-4110KVBK
● Cut to fit 15 to 28 mm	28 mm	LPSCA28-4115KV
● Pre-cut to 20 mm	28 mm	LPSCA28-4120KV
● Pre-cut to 25 mm	28 mm	LPSCA28-4125KV
● Cut to fit 15 to 40 mm	40 mm	LPSCA40-4115KVBK
● Cut to fit 15 to 40 mm	40 mm	LPSCA40-4115KV
● Pre-cut to 30 mm	40 mm	LPSCA40-4130KV
● Pre-cut to 35 mm	40 mm	LPSCA40-4135KV
● Cut to fit 15 to 47 mm	47 mm	LPSCA47-4115KVBK
● Cut to fit 15 to 47 mm	47 mm	LPSCA47-4115KV
● Pre-cut to 40 mm	47 mm	LPSCA47-4140KV
● Pre-cut to 45 mm	47 mm	LPSCA47-4145KV

● Black cover ● Beige cover



Alginate



Alginate

Drainable Standard Alginate Soft Convex

7 mm plateau height

SCA25-4110KVBK - SCA45-4140KV Window, opaque, filter, soft cover both sides, soft touch closure.

Pack Size: 20

Hole size	Plateau	Window
● Cut to fit 10 to 24 mm	25 mm	SCA25-4110KVBK
● Cut to fit 20 to 24 mm	25 mm	SCA25-4120KV
● Cut to fit 15 to 34 mm	35 mm	SCA35-4115KVBK
● Cut to fit 20 to 34 mm	35 mm	SCA35-4120KV
● Pre-cut to 25 mm	35 mm	SCA35-4125KV
● Pre-cut to 30 mm	35 mm	SCA35-4130KV
● Cut to fit 15 to 44 mm	45 mm	SCA45-4115KVBK
● Cut to fit 20 to 44 mm	45 mm	SCA45-4120KV
● Pre-cut to 35 mm	45 mm	SCA45-4135KV
● Pre-cut to 40 mm	45 mm	SCA45-4140KV

● Black cover ● Beige cover



Primo hydrocolloid

Drainable Standard Primo Hydrocolloid

HP10D20 - HP10D60 Window, opaque, filter, soft cover both sides, soft touch closure.

HP10D20C Clear, filter, soft cover one side, soft touch closure.

Pack Size: 30

Hole size	Window	Clear
Cut to fit 10 to 50 mm	HP10D20	HP10D20C
Pre-cut to 25 mm	HP10D25	
Pre-cut to 30 mm	HP10D30	
Pre-cut to 35 mm	HP10D35	
Pre-cut to 40 mm	HP10D40	
Pre-cut to 45 mm	HP10D45	
Pre-cut to 50 mm	HP10D50	
Pre-cut to 55 mm	HP10D55	
Pre-cut to 60 mm	HP10D60	



Primo hydrocolloid

Drainable Midi Primo Hydrocolloid

HP8D20 - HP8D60 Window, opaque, filter, soft cover both sides, soft touch closure.

HP8D20C Clear, filter, soft cover one side, soft touch closure.

Pack Size: 30

Hole size	Window	Clear
Cut to fit 20 to 60 mm	HP8D20	HP8D20C
Pre-cut to 25 mm	HP8D25	
Pre-cut to 30 mm	HP8D30	
Pre-cut to 35 mm	HP8D35	
Pre-cut to 40 mm	HP8D40	
Pre-cut to 45 mm	HP8D45	
Pre-cut to 50 mm	HP8D50	
Pre-cut to 55 mm	HP8D55	
Pre-cut to 60 mm	HP8D60	



Primo hydrocolloid

Drainable Standard Primo Hydrocolloid High Output

HVP20-60 Window, opaque, filter, soft cover both sides, wide bung closure.

Pack Size: 10

Hole size	Window
Cut to fit 20 to 60 mm	HVP20-60



Drainable Standard Primo Hydrocolloid High Output Soft Convex

7 mm plateau height

SCHV25-4115K - SCHV45-4115K Window, opaque, filter, soft cover both sides, wide bung closure.
Pack Size: 10

Hole size	Plateau	Window
Cut to fit 15 to 24 mm	25 mm	SCHV25-4115K
Cut to fit 15 to 34 mm	35 mm	SCHV35-4115K
Cut to fit 15 to 44 mm	45 mm	SCHV45-4115K



Primo hydrocolloid



Drainable Midi Primo Hydrocolloid Low Profile Soft Convex

3.5 mm plateau height

LPSC28-4315KV - LPSC47-4345KV Window, opaque, filter, soft cover both sides, soft touch closure.
LPSC28-4315KVC - LPSC47-4315KVC Clear, filter, soft cover one side, soft touch closure.
Pack Size: 20

Hole size	Plateau	Window	Clear
Cut to fit 15 to 28 mm	28 mm	LPSC28-4315KV	LPSC28-4315KVC
Pre-cut to 20 mm	28 mm	LPSC28-4320KV	
Pre-cut to 25 mm	28 mm	LPSC28-4325KV	
Cut to fit 15 to 40 mm	40 mm	LPSC40-4315KV	LPSC40-4315KVC
Pre-cut to 30 mm	40 mm	LPSC40-4330KV	
Pre-cut to 35 mm	40 mm	LPSC40-4335KV	
Cut to fit 15 to 47 mm	47 mm	LPSC47-4315KV	LPSC47-4315KVC
Pre-cut to 40 mm	47 mm	LPSC47-4340KV	
Pre-cut to 45 mm	47 mm	LPSC47-4345KV	



Primo hydrocolloid



Drainable Standard Primo Hydrocolloid Low Profile Soft Convex

3.5 mm plateau height

LPSC28-4110KVBK - LPSC47-4145KV Window, opaque, filter, soft cover both sides, soft touch closure.
LPSC28-4115KVC - LPSC47-4115KVC Clear, filter, soft cover one side, soft touch closure.
Pack Size: 20

Hole size	Plateau	Window	Clear
● Cut to fit 10 to 28 mm	28 mm	LPSC28-4110KVBK	
● Cut to fit 15 to 28 mm	28 mm	LPSC28-4115KV	LPSC28-4115KVC
● Pre-cut to 20 mm	28 mm	LPSC28-4120KV	
● Pre-cut to 25 mm	28 mm	LPSC28-4125KV	
● Cut to fit 15 to 40 mm	40 mm	LPSC40-4115KVBK	
● Cut to fit 15 to 40 mm	40 mm	LPSC40-4115KV	LPSC40-4115KVC
● Pre-cut to 30 mm	40 mm	LPSC40-4130KV	
● Pre-cut to 35 mm	40 mm	LPSC40-4135KV	
● Cut to fit 15 to 47 mm	47 mm	LPSC47-4115KVBK	
● Cut to fit 15 to 47 mm	47 mm	LPSC47-4115KV	LPSC47-4115KVC
● Pre-cut to 40 mm	47 mm	LPSC47-4140KV	
● Pre-cut to 45 mm	47 mm	LPSC47-4145KV	

● Black cover ● Beige cover



Primo hydrocolloid



Drainable Midi Primo Hydrocolloid Soft Convex

7 mm plateau height


SC25-4315KV - SC45-4340KV Window, opaque, filter, soft cover both sides, soft touch closure.
 SC25-4315KVC - SC45-4315KVC Clear, filter, soft cover one side, soft touch closure.
 Pack Size: 20



Primo hydrocolloid

Hole size	Plateau	Window
Cut to fit 15 to 24 mm	25 mm	SC25-4315KV
Pre-cut to 20 mm	25 mm	SC25-4320KV
Cut to fit 15 to 34 mm	35 mm	SC35-4315KV
Pre-cut to 25 mm	35 mm	SC35-4325KV
Pre-cut to 30 mm	35 mm	SC35-4330KV
Cut to fit 15 to 44 mm	45 mm	SC45-4315KV
Pre-cut to 35 mm	45 mm	SC45-4335KV
Pre-cut to 40 mm	45 mm	SC45-4340KV

Hole size	Plateau	Clear
Cut to fit 15 to 24 mm	25 mm	SC25-4315KVC
Cut to fit 15 to 34 mm	35 mm	SC35-4315KVC
Cut to fit 15 to 44 mm	45 mm	SC45-4315KVC



Drainable Standard Primo Hydrocolloid Soft Convex

7 mm plateau height

C25-4110KVBK - SC57-4120KV Window, opaque, filter, soft cover both sides, soft touch closure.
 Pack Size: 20



Primo hydrocolloid

Hole size	Plateau	Window
● Cut to fit 10 to 24 mm	25 mm	SC25-4110KVBK
● Cut to fit 20 to 24 mm	25 mm	SC25-4120KV
● Cut to fit 15 to 34 mm	35 mm	SC35-4115KVBK
● Cut to fit 20 to 34 mm	35 mm	SC35-4120KV
● Pre-cut to 25 mm	35 mm	SC35-4125KV
● Pre-cut to 30 mm	35 mm	SC35-4130KV
● Cut to fit 15 to 44 mm	45 mm	SC45-4115KVBK
● Cut to fit 20 to 44 mm	45 mm	SC45-4120KV
● Pre-cut to 35 mm	45 mm	SC45-4135KV
● Pre-cut to 40 mm	45 mm	SC45-4140KV

Pack Size: 10

Hole size	Plateau	Window
● Cut to fit 20 to 57 mm	57 mm	SC57-4120KV

● Black cover ● Beige cover



Drainable Standard Microskin

MV-4120K Window, filter, soft cover both sides, soft touch closure.
 Pack Size: 30



Microskin

Hole size	Standard
Cut to fit 20 to 50 mm	MV-4120K

Urostomy

Continuous drainage requires solutions that provide dependable flow while protecting the skin around the stoma. Comfort, security and simplicity are key in everyday urostomy care.

Unora urostomy solutions combine reliable drainage with flexible skin barriers and adaptability for drainage systems, helping protect the skin and maintain confidence.

Our extended urostomy range is shaped for maximum capacity with minimum bulk. The unique non-twisting pouch design allows for easier emptying. For the outlet, there is a flexible choice between a soft bung or tap closure for comfortable wear in everyday life.



Primo hydrocolloid

Dual-layer hydrocolloid technology engineered to enhance adhesion and security while remaining flexible and comfortable.



Microskin

Ultra-thin, flexible skin protection that adapts closely to the body and supports natural movement.

Product overview



Primo hydrocolloid

Urostomy Standard Primo Hydrocolloid	46
Urostomy Standard Primo Hydrocolloid Low Profile Soft Convex	48
Urostomy Standard Primo Hydrocolloid Soft Convex	50



Microskin

Urostomy Standard Microskin	52
-----------------------------	----

Starter pack



Clinic/ Community

Our closed bag starter packs are perfect for supporting patients in finding the bag that suits them best. Available in all standard and maxi cut to fit closed bags, each pack contains 10 bags of the selected product.

Full instructions for use are included to provide guidance on application and removal, along with a sizing guide to ensure the bag wafer is cut to the correct size. The cover card details all the starter packs available and may be used as a reference tool to choose an alternative bag, if the initial one chosen does not suit the patient.

Pack contains

- 5 x Bags of chosen product
- 1 x Night drainage bag
- 2 x Connectors
- 1 x Belt
- 1 x IFU
- 1 x Product code information sheet

Products that are available to order as a starter pack are indicated in the product details on the following pages.

Please contact our customer service team on **0800 592 786** or email samples@unora.co.uk to place an order.



Urostomy Standard Primo Hydrocolloid

GC1500 Window, opaque, soft cover both sides, soft bung closure, capacity 420ml.
Pack Size: 20

Hole size	Window
● Cut to fit 15 – 60 mm	GC1500

- Available as a Clinic/Community pack. Call us to order.



Primo hydrocolloid



Urostomy Standard Primo Hydrocolloid


GC1500C Clear, soft cover one side, soft bung closure, capacity 420ml.
Pack Size: 20

Hole size	Clear
● Cut to fit 15 – 60 mm	GC1500C

- Available as a Clinic/Community pack. Call us to order.



Primo hydrocolloid



Urostomy Standard Primo Hydrocolloid


GC1500T Window, opaque, soft cover both sides, tap closure, capacity 420ml.
Pack Size: 20

Hole size	Window
● Cut to fit 15 – 60 mm	GC1500T

- Available as a Clinic/Community pack. Call us to order.



Primo hydrocolloid



Urostomy Standard Primo Hydrocolloid

GC1500CT Clear, soft cover both sides, tap closure, capacity 420ml.
Pack Size: 20

Hole size	Clear
● Cut to fit 15 – 60 mm	GC1500CT

- Available as a Clinic/Community pack. Call us to order.



Primo hydrocolloid



Urostomy Standard Primo Hydrocolloid Low Profile Soft Convex

3.5 mm
plateau height

LPSCU28-0015 - LPSCU47-0045 Window, opaque, soft cover both sides, bung closure, capacity 420ml. Pack Size: 20

Hole size	Plateau	Window
● Cut to fit 15 to 28 mm	28 mm	LPSCU28-0015
Pre-cut to 20 mm	28 mm	LPSCU28-0020
Pre-cut to 25 mm	28 mm	LPSCU28-0025
● Cut to fit 15 to 40 mm	40 mm	LPSCU40-0015
Pre-cut to 30 mm	40 mm	LPSCU40-0030
Pre-cut to 35 mm	40 mm	LPSCU40-0035
● Cut to fit 15 to 47 mm	47 mm	LPSCU47-0015
Pre-cut to 40 mm	47 mm	LPSCU47-0040
Pre-cut to 45 mm	47 mm	LPSCU47-0045



Primo hydrocolloid

● Available as a Clinic/Community pack.
Call us to order.



Urostomy Standard Primo Hydrocolloid Low Profile Soft Convex

3.5 mm
plateau height

LPSCU28-0015T - LPSCU47-0045T Window, opaque, soft cover both sides, tap closure, capacity 420ml. Pack Size: 20

Hole size	Plateau	Window
● Cut to fit 15 to 28 mm	28 mm	LPSCU28-0015T
Pre-cut to 20 mm	28 mm	LPSCU28-0020T
Pre-cut to 25 mm	28 mm	LPSCU28-0025T
● Cut to fit 15 to 40 mm	40 mm	LPSCU40-0015T
Pre-cut to 30 mm	40 mm	LPSCU40-0030T
Pre-cut to 35 mm	40 mm	LPSCU40-0035T
● Cut to fit 15 to 47 mm	47 mm	LPSCU47-0015T
Pre-cut to 40 mm	47 mm	LPSCU47-0040T
Pre-cut to 45 mm	47 mm	LPSCU47-0045T



Primo hydrocolloid

● Available as a Clinic/Community pack.
Call us to order.



Urostomy Standard Primo Hydrocolloid Low Profile Soft Convex

3.5 mm
plateau height


LPSCU28-0015C - LPSCU47-0015C Clear, soft cover one side, soft bung closure, capacity 420ml. Pack Size: 20

Hole size	Plateau	Clear
● Cut to fit 15 to 28 mm	28 mm	LPSCU28-0015C
● Cut to fit 15 to 40 mm	40 mm	LPSCU40-0015C
● Cut to fit 15 to 47 mm	47 mm	LPSCU47-0015C



Primo hydrocolloid

● Available as a Clinic/Community pack.
Call us to order.



Urostomy Standard Primo Hydrocolloid Low Profile Soft Convex

3.5 mm
plateau height

LPSCU28-0015CT - LPSCU47-0015CT Clear, soft cover one side, tap closure, capacity 420ml. Pack Size: 20

Hole size	Plateau	Clear
● Cut to fit 15 to 28 mm	28 mm	LPSCU28-0015CT
● Cut to fit 15 to 40 mm	40 mm	LPSCU40-0015CT
● Cut to fit 15 to 47 mm	47 mm	LPSCU47-0015CT



Primo hydrocolloid

● Available as a Clinic/Community pack.
Call us to order.



Urostomy Standard Primo Hydrocolloid Soft Convex

7 mm
plateau height


SCU25-0015 - SCU45-0040 Window, opaque,
soft cover both sides, soft bung closure, capacity 420ml.
Pack Size: 20

Hole size	Plateau	Window
● Cut to fit 15 to 24 mm	25 mm	SCU25-0015
Pre-cut to 20 mm	25 mm	SCU25-0020
● Cut to fit 15 to 34 mm	35 mm	SCU35-0015
Pre-cut to 25 mm	35 mm	SCU35-0025
Pre-cut to 30 mm	35 mm	SCU35-0030
● Cut to fit 15 to 44 mm	45 mm	SCU45-0015
Pre-cut to 35 mm	45 mm	SCU45-0035
Pre-cut to 40 mm	45 mm	SCU45-0040



Primo hydrocolloid

- Available as a Clinic/Community pack.
Call us to order.



Urostomy Standard Primo Hydrocolloid Soft Convex

7 mm
plateau height

SCU25-0015C - SCU45-0015C Clear, soft cover one side,
soft bung closure, capacity 420ml.
Pack Size: 20

Hole size	Plateau	Clear
● Cut to fit 15 to 24 mm	25 mm	SCU25-0015C
● Cut to fit 15 to 34 mm	35 mm	SCU35-0015C
● Cut to fit 15 to 44 mm	45 mm	SCU45-0015C



Primo hydrocolloid

- Available as a Clinic/Community pack.
Call us to order.



Urostomy Standard Primo Hydrocolloid Soft Convex

7 mm
plateau height


SCU25-0015T - SCU45-0040T Window, opaque,
soft cover both sides, tap closure, capacity 420ml.
Pack Size: 20

Hole size	Plateau	Window
● Cut to fit 15 to 24 mm	25 mm	SCU25-0015T
Pre-cut to 20 mm	25 mm	SCU25-0020T
● Cut to fit 15 to 34 mm	35 mm	SCU35-0015T
Pre-cut to 25 mm	35 mm	SCU35-0025T
Pre-cut to 30 mm	35 mm	SCU35-0030T
● Cut to fit 15 to 44 mm	45 mm	SCU45-0015T
Pre-cut to 35 mm	45 mm	SCU45-0035T
Pre-cut to 40 mm	45 mm	SCU45-0040T



Primo hydrocolloid

- Available as a Clinic/Community pack.
Call us to order.



Urostomy Standard Primo Hydrocolloid Soft Convex

7 mm
plateau height

SCU25-0015CT - SCU45-0015CT Clear, soft cover one side,
tap closure, capacity 420ml.
Pack Size: 20

Hole size	Plateau	Clear
● Cut to fit 15 to 24 mm	25 mm	SCU25-0015CT
● Cut to fit 15 to 34 mm	35 mm	SCU35-0015CT
● Cut to fit 15 to 44 mm	45 mm	SCU45-0015CT



Primo hydrocolloid

- Available as a Clinic/Community pack.
Call us to order.



Microskin

Urostomy Standard Microskin

MJ1300T - MJ3800T Window, opaque, soft cover both sides, tap closure.

Pack Size: 20

Hole size	Window
Cut to fit 10 to 50 mm	MJ1300T
Pre-cut to 16 mm	MJ1600T
Pre-cut to 19 mm	MJ1900T
Pre-cut to 22 mm	MJ2200T
Pre-cut to 25 mm	MJ2500T
Pre-cut to 32 mm	MJ3200T
Pre-cut to 38 mm	MJ3800T



Microskin

Urostomy Standard Microskin

MJ1300 - MJ3800 Window, opaque, soft cover both sides, soft bung closure.

Pack Size: 20

Hole size	Window
Cut to fit 10 to 50 mm	MJ1300
Pre-cut to 16 mm	MJ1600
Pre-cut to 19 mm	MJ1900
Pre-cut to 22 mm	MJ2200
Pre-cut to 25 mm	MJ2500
Pre-cut to 32 mm	MJ3200
Pre-cut to 38 mm	MJ3800



Wound management

Some situations call for solutions beyond everyday pouching systems. From neonatal care to complex drainage needs, specialised products help support safe and effective stoma management when standard options are not enough.

Our special supply range has been developed for medical professionals and offers specialised solutions that address these specific challenges, helping protect delicate skin, support clinical care and provide reassurance when it matters most.

All our wound managers come with the option of a port, enabling easy access to the wound.

Option Manager	56
Option Wound Manager Extra Small	57
Option Wound Manager Small	57
Option Wound Manager Oblong	57
Option Wound Manager Medium	58
Option Wound Manager Square	58
Option Wound Manager Large	59
Option Wound Manager Extra Large	59

Option Manager

OM-S Clear, soft cover one side, cap closure.

Approx capacity 50ml.

OM-M Clear, soft cover one side, cap closure.

Approx capacity 150ml.

OM-L Clear, soft cover one side, cap closure.

Approx capacity 300ml.

Pack Size: 30



Hydrocolloid

Hole size	Window
Cuts up to 30 mm	OM-S
Cuts up to 50 mm	OM-M
Cuts up to 50 mm	OM-L



Hydrocolloid

Option Wound Manager Extra Small

WWM-XS - WWM-XSA Hole 10 mm, film: clear.

Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 90 mm wide x 180 mm high	WWM-XS	WWM-XSA

(Oval shape 90 mm x 180 mm at maximum point)



Hydrocolloid

Option Wound Manager Small

WWM-S - WWM-SA Hole 10 mm, film: clear.

Pack Size: 10

Hole Size	Standard	Access Port
Suitable for horizontal wounds up to 240 mm wide x 155 mm high	WWM-S	WWM-SA



Hydrocolloid

Option Wound Manager Oblong

WWM-O - WWM-OA Hole 10 mm, film: clear.

Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 60 mm wide x 190 mm high	WWM-O	WWM-OA



Option Wound Manager Medium

WWM-M - WWM-MA Hole 10 mm, film: clear.

Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 85 mm wide x 250 mm high	WWM-M	WWM-MA



Hydrocolloid



Option Wound Manager Large

WWM-L - WWM-LA Hole 10 mm, film: clear.

Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 160 mm wide x 260 mm high	WWM-L	WWM-LA



Hydrocolloid



Option Wound Manager Square

WWM-SQ - WWM-SQA Hole 10 mm, film: clear.

Pack Size: 10

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 160 mm wide x 200 mm high	WWM-SQ	WWM-SQA



Hydrocolloid



Option Wound Manager Extra Large*

WWM-XL - WWM-XLA Hole 10 mm, film: clear.

Pack Size: 5

Hole Size	Standard	Access Port
Suitable for vertical wounds up to 230 mm wide x 380 mm high	WWM-XL	WWM-XLA



Hydrocolloid

* Option Wound Manager Extra Large is not available on prescription.

Accessories

Good stoma care extends beyond the pouching system itself. Small details – protecting the skin, improving adhesion or supporting comfortable wear – can make a meaningful difference in everyday life.

Our accessories bring together practical solutions designed to complement pouching systems and support peristomal skin health, helping create a reliable and comfortable foundation for daily stoma care.

EasyLease Sting Free Adhesive Remover Spray	62
Elite® Sting Free Adhesive Remover	62
Elisse® Sting Free Skin Barrier	62
Option Seals	63
Genii™ Ostomy Seals	63
Genii™ Convex Ostomy Seals	63
Curvies	64
Curvies Plus	64
Curvies Hydro	64
Genii™ Silicone Flange Extenders	64
Pearls® Gelling & Odour Control Sachets	65
Ostomy Belt	65
Urostomy Drainage Adaptor	65
Hydrocolloid Stoma Cap with Soft Cover and Filter	65



EasyLease Sting Free Adhesive Remover Spray

No sting medical adhesive remover with antiseptic tea tree oil.

	Code
50ml Adhesive Remover Aerosol	1408



Option Seals

High tack and absorbent.
Pack of 10

	Code
Slim Seals - Slim and stretchy, ideal for skin protection.	WM-500
Seals - Slightly thicker than the slim seals for added depth.	WM-100



Elite® Sting Free Adhesive Remover

100% silicone formulation.
Hypoallergenic and sting free spray, non-irritant adhesive remover.
Also available in a wipe.

	Code
50ml spray can	TR201
Pack of 30 wipes	TR202



Genii™ Ostomy Seals

Protects the skin immediately around the stoma from irritation and breakdown
Pack Size: 30

	Code
Small - 20 mm-35 mm	TR1020
Medium - 35 mm-50 mm	TR1028
Large - 50 mm-80 mm	TR1035



Elisse® Sting Free Skin Barrier

Hypoallergenic silicone formula that builds a protective barrier between the ostomy wafer or base plate and the skin.
Also available in a wipe.

	Code
50ml spray can	TR204
Pack of 30 wipes	TR203



Genii™ Convex Ostomy Seals

Soft and mildly convex to provide the best fit to the body shape. Can be applied convex side up or down.
Pack Size: 10

	Code
Small - 20 mm-35 mm	TR1020cx
Large - 30 mm-50 mm	TR1030cx



Curvies

Secure and flexible silicone adhesive flange extenders. 110 mm max width. CUR001 Pack Size: 30

	Code
Curvies	CUR001



Curvies Plus

Secure and flexible silicone adhesive, large flange extenders. 155 mm max width. CUR005 Pack Size: 20

	Code
Curvies Plus	CUR005



Curvies Hydro

Kind to your skin, Hydrocolloid adhesive, large flange extenders. 155 mm max width. HYD007 Pack Size: 20

	Code
Curvies Hydro	HYD007



Genii™ Silicone Flange Extenders

Provides a soft, flexible, invisible border which conforms to the body. TR1060 Pack Size: 20

	Code
Genii™ Flange Extenders	TR1060



Pearls® Gelling & Odour Control Sachets

Super absorbent gelling crystals to be placed inside the ostomy bag.

	Code
Tub of 100 sachets	TR205




Ostomy Belt

	Code
1 per pack	OCB02



Urostomy Drainage Adaptor

	Code
5 per pack	UDA



Hydrocolloid Stoma Cap with Soft Cover and Filter

1420K Filter, opaque, soft cover both sides. Pack Size: 50

Hole Size	Soft Cover - Opaque
Cut to fit 20 to 45 mm	1420K

Here for what matters

Unora Home Delivery provides patients with a reliable and compassionate home delivery service for all stoma prescription products.

For free help and advice, call us on 0800 592 786 or email us at support@unorahomedelivery.co.uk

If you are a healthcare professional registering a new patient please email: priority@unorahomedelivery.co.uk



Unora Home Delivery®

Unora Healthcare Ltd
Unit B
Barn Way
Lodge Farm Industrial Estate
Northampton
NN5 7UW

contact@unora.co.uk
www.unora.co.uk

unora